1. Intended use / Indication

otosil® ES-1.1 from pro3dure is a 1:1 Silicone material for manufacturing of earmolds and hearing protections.

2. Contraindication

otosil® ES-1.1 is contraindicated ...
1. ... If a patient is known to be allergic to any of the ingredients.
2. ... for every application that is not part of the indication (see above).

3. Patient Target Group

Persons for whom an earmold is to be created.

4. Intended Users

Earmold laboratory, ENT specialist, hearing aid acoustician

5. Requirements

Use by specialists only.

6. Material

Contains: Functionalized polydimethylsiloxanes, crosslinkers, modifiers and additives.

7. Manufacturing process (fig. 1-6)

Mixing and dosing

1. Place the cartridge into the dispensing gun.
2. Unscrew cap. In order to ensure a proper flow and mixing ratio from both orifices extrude carefully a small amount of the material.
3. Insert the static mixer into the guided grooves of cartridge. Rotate the mixing cannula anti-clockwise to a stop. If necessary, mount an Intra-tip. For reducing forces during dispensing, static mixers of a larger diameter (9 mm) can be used alternatively. Dispenser is now ready to mix the 2 component earmold silicone in any amount needed.

Application (PnP-method)

4. With the otosil® ES-1.1 material all commercially available materials for negative forms like plasters, gels or 3D-printing resins can be used for the manufacturing of the negative form.
5. If necessary, coat the negative form with a separation liquid (e.g. alginate based for plaster). 3D-printed cast forms have to be cleaned carefully before use in order to avoid inhibition of the silicone reaction. Inject the otosil® ES-1.1 material slowly and bubblefree into the generated negative form. In case of narrow structures (e.g. auditory canals), we recommend using a thin mixing cannula. Put the negative mold with the injected material into a pressure pot for about 60 min at 122 °F (50 °C).

Final surface treatment and lacquering

6. The silicone earmold can now be shaped by using special cutters and grinding sleeves. Roughen the surface with corundum paper (grit 180) before varnishing. Clean the surface and coat it with a silicone lacquer according to the special instructions of the lacquer.

8. Important working hints


In case of the generative fabricated castforms a careful preparation as well as a thorough cleaning will ensure better and more reliable results. Ensure that all traces of residual resin is removed from the negative form. The cured earmold material otosil® ES-1.1 is chemically inert. Avoid spots on clothing. We recommend standard gloves made of nitrile or polyethylene.

9. Notice

Undesired reactions of the immune system like allergies, irritations can't be absolutely excluded. In case of doubt, we recommend to contact your consulting physician and make an allergy test before the application of the material. According to the EU Medical Devices Regulation, users/patients are obliged to report serious events with a medical device to the manufacturer and to the competent authority of the country in which they occurred. Dispose of contents/container according to official regulations.

**Earmold silicone
ES-1.1 40 Shore**

Instructions for use

Physical properties*:

otosil® ES-1.1 40 Shore

- Setting time (50 °C) / min.: ≤ 15**
- Processing time (23 °C) / min.: ≥ 4**
- Tensile strength / MPa (ISO 37:2017): ≥ 2**
- Elongation at break / % (ISO 37:2017): ≥ 400 **
- Shore hardness A (DIN ISO 48-4:2021): 40±5**

Ordering information:

otosil® ES-1.1 40 Shore

Standard package:

Cartridge
8 x 48 ml
(à 2 x 24 ml A + B)

clear-transparent
REF: A0032200

reddish-transparent
REF: A0031603

dark blue-opaque
REF: A0031828

red-opaque
REF: A0032027

Additional colors on request.

Cannula:
A System
5,4 mm
transparent
50 pieces

REF: A4020142

Dispenser:
DMA50

REF: A4020128

* These data come from measurements of a representative sample that was determined as part of our quality assurance.
** According to internal design and requirements specifications.

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pro3dure medical GmbH
pro3dure medical LLC