### 1. Intended use / Indication

The **printo<sup>®</sup> GR-10 guide | MSI** resin is a light-curable polymerizable resin intended to be used in conjunction with extra-oral curing light equipment for the fabrication, by additive manufacturing, of orthodontic and dental objects such as splints.

### 2. Contraindication

**printo<sup>®</sup> GR-10 guide | MSI** is contraindicated ...

1. ... if a patient is known to be allergic to any of the ingredients.
2. ... for bruxism splints and splints with wall thickness ≤ 1.5mm.
3. ... for denture bases.
4. ... for every application that is not part of the indication (see above).

### 3. Product Description

The **printo<sup>®</sup> GR-10 guide | MSI** resin includes a lactam-based, biomimetic solution for multi-species inhibition\*. This reduces the formation of undesirable biofilms on orthodontic appliances without a microcidal effect. \*c. albicans / s. mutans

### 4. Patient Target Group

Persons being treated in the context of a dental procedure.

### 5. Intended Users

Dentist, dental technicians

### 6. Requirements

**Software – information available from:**

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**Hardware (3D Printing) – information available from:**

ASIGA · 2/19-21 Bourke Road, Alexandria, NSW 2015, Australia

**Hardware (Post Curing) – information available from:**

pro3dure medical (CD-1, CD-2) · Am Burgberg 13, D-58642 Iserlohn, Germany

Additional specifications on [www.pro3dure.com](http://www.pro3dure.com).

### 7. Material

**printo<sup>®</sup> GR-10 guide | MSI** consists of functional (meth)acrylic resins, initiators, dyes and stabilisers.

### 8. Geometric presetting

Minimum wall thickness: 1.5 mm

### 9. Material parameters

Radiation penetration depth controllable by exposure time

50 µm

100 µm

### 10. Manufacturing process (fig. 1–10)

1. Prepare data (CAD & build preparation).
2. Choose process parameters (build style, etc.).
3. Transfer prepared data to 3D printing device.
4. Prepare 3D printing – shake bottle.
5. Fill resin tank of 3D printing device.
6. Build the parts.
7. Clean parts (with IPA ≥ 97 % or equivalent cleaning agent) approx. 4 min. in an ultrasonic bath or equivalent device – precleaning recommended.
8. Dry parts, until there are no residues of IPA or equivalent cleaning agent.
9. Post curing (4 min.): inert atmosphere recommended (use adequate light curing device).
10. Finish parts.

### 11. Finishing processes

Polishing


### 12. Notice

Please follow the instructions provided by the software manufacturer in terms of parameter settings and design recommendations. Please follow the instructions provided by the hardware manufacturer in terms of parameter settings/printing & post curing recommendations. To avoid detrimental effects on material quality do not expose the liquid material to irradiation under any circumstances. Deviations from the described manufacturing processes or storing conditions may lead to different mechanical and optical properties of the material. According to the EU Medical Devices Regulation, users/patients are obliged to report serious events with a medical device to the manufacturer and to the competent authority of the country in which they occurred. **Caution:** The lot number and the best before date are indicated on each packaging. In case of claims please always indicate the lot number of the product. Do not use the product after expiry of the best before date. Avoid breathing dust/fume/gas/mist/vapours/spray. Wash hands thoroughly after handling. Use only outdoors or in a well-ventilated area. Contaminated work clothing should not be allowed out of the workplace. Wear protective gloves/protective clothing/eye protection/face protection. IF ON SKIN: Wash with plenty of soap and water. IF INHALED: Remove person to fresh air and keep comfortable for breathing. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Call a POISON CENTER or doctor/physician if you feel unwell. Specific treatment (see information on this label). If skin irritation or rash occurs: Get medical advice/attention. If eye irritation persists: Get medical advice/attention. Take off contaminated clothing. Store in a well-ventilated place. Keep container tightly closed. Store locked up. Dispose of contents/container according to official regulations.

### 13. Hazard statements

Causes skin irritation. May cause an allergic skin reaction. Causes serious eye irritation. May cause respiratory irritation.

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
**Distributor (US):**  
pro3dure medical LLC  
9825 Valley View Road · Eden Prairie, MN 55344  
Phone: +1 (952)-426-1928

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