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THERMEO® SO Fräsrohlinge dienen zur Herstellung funktioneller Ästhetikschienen (Snap-On-Schienen).

#### 2. Produktbeschreibung

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
Distributor (US):  
pro3dure medical LLC  
9825 Valley View Road · Eden Prairie, MN 55344  
Phone: +1 (952)-426-1928

info@pro3dure.com · www.pro3dure.com · Made in Germany

### Instruction for use

- GB 1. Usage
- US THERMEO® SO milling blank is used for the production of dental snap on splints.

### 2. Product description

THERMEO® SO milling blanks are used to manufacture dental snap on splints with thermomemory effect. It is characterized by its thermoplastic flexibility, which enables a highly precised adaptation to the tooth situation. It offers stress-free wearing comfort for the patients. On top of that, the self-adjusting THERMEO® SO splints are extremely resilient to breakage and present an optimized translucency within the tooth colors A1, A2 and A3. The industrial production process of THERMEO® SO provides an exceptional material homogeneity and thereby an outstanding long-term stability. The usage of THERMEO® SO milling blank with help of CAD/CAM technology ensures a safe production process, because the mixing failure (e.g. through hand mixing) is not applied. Moreover, the smell, which is caused by the mixing process, can be avoided. Indications for THERMEO® SO splints: long-term temporaries for concealing aesthetic deficiencies, protecting and releasing splints, splints for diagnostical bite elevation.

### 3. Processing

- THERMEO® SO milling blanks should be only used in the CAD/CAM production of dental snap on splints and can be used in all common CAD/CAM systems (Fig. 1). Using the material for other indications is contraindicated. There should be only hard metal cutters used during the processing.
- THERMEO® SO snap on splints are designed and milled using CAD/CAM technology by qualified persons (dental technicians). After the finishing, the occlusal thickness of 0.2 mm needs to be reached. Any corners and edges must be rounded off. To remove the THERMEO® SO snap on splints from the blanks, a finely toothed hard metal cutter or a suitable cutting tool is recommended to be used. Carefully cut through the retaining strips without pressure.
- In order to avoid the plaque accumulations subsequent polishing of the splint is essential. Please prepolish with suitable polishes and a goat hair brush. The high-luster polishing should be done with a corresponding polishing compound (Fig. 2). As far as possible, please avoid heat generation during the polishing and finishing of the THERMEO® SO splint. This guarantees the optimal fitting of the splint.
- Please note:** Upon completion, use an ultrasonic bath with water at a maximum temperature of approx. 35 °C for cleaning. (Important: steam cleaning devices are not recommended, because the product at higher temperature and pressure can be deformed.) The processing recommendations of the manufacturer must be taken in consideration.


### 4. Milling parameters THERMEO® SO

Roughing:	Finishing:	Rest roughing:
■ 2 mm ball radius cutter – single cutter	■ 2 mm ball radius cutter – single cutter	■ 1 mm ball radius cutter – single cutter
■ Speed: 22.000 rpm	■ Speed: 22.000 rpm	■ Speed: 28.000 rpm
■ Radial path distance: 0,6 mm	■ Radial path distance: 0,1 mm	■ Radial path distance: 0,1 mm
■ Z depth increase: 0,6 mm	■ Z depth increase: 0,1 mm	■ Z depth increase: 0,1 mm
■ Feed: 1.600 mm/min	■ Feed: 1.600 mm/min	■ Feed: 1.200 mm/min

The irrigation nozzle on the spindle should be aligned with the tool tip. In this occasion it should be noted that a uni-form programming length is respected or the length of the 2 mm tool is estimated. The cooling medium air is sufficient.

### 5. Notes


- Incorrect use of milling tools, polishing brushes, steam cleaner and water baths can overheat the material exposed and thus the properties are adversely affected. Deviation from the described manufacturing process may lead to different mechanical and optical properties of the material.
- The lot number and the best before date are indicated on each THERMEO® SO packaging. In case of claims please always indicate the lot number of the product.
- Do not use the product after expiration date.

Milling blank for the production of snap on splints

#### Technical data:

- Colour: A1, A2, A3
- Density: ca. 1.1 g/ml
- Flexural strength: 23 °C (MPa) > 20
- Flexural strength: 37 °C (MPa) < 20
- Storage:

#### Ordering information:

THERMEO® SO-blank  
d = 20 mm  
∅ = 98,5 mm

- A1: item no.: D1010125
- A2: item no.: D1010126
- A3: item no.: D1010127


These data result from measurements of a representative sample, which were determined within the scope of our quality assurance.

- The use of the material for the production of conventional splints is contraindicated, as is the production of denture bases. Undesirable biological reactions (e.g. allergies to material components) may occur in individual cases. Clarify the use of THERMEO® SO splints with your dentist or orthodontist in case of a known allergic response to components of THERMEO® SO.
- Clean THERMEO® SO splints under cold running water with a toothbrush. Disinfect the splint with non-alcoholic detergents before first use.
- According to the EU Medical Devices Regulation, users/patients are obliged to report serious events with a medical device to the manufacturer and to the competent authority of the country in which they occurred.

### 6. Safety advice

- Restorations are custom-made devices in accordance with the Directive 93/42/EEC. Record the BATCH no. for every procedure requiring identification of the material.
- **Storage:** No special measures are required.
- **Disposal:** Dental objects made of THERMEO® SO are not water soluble, are inert, and do not present a hazard for ground water. As such, they can be disposed of as domestic waste without any special precautions if not stated otherwise in your local disposal guidelines.
- **Warranty:** Our recommendations concerning the application technique, regardless of whether they are communicated in writing, orally, or by means of practical instructions, are based on our own experiences and tests. As such, they are intended as guidelines only. We are continually striving to improve our products. Consequently, we reserve the right to make changes to their design and composition.
- **Hazard Information:** Processing of THERMEO® SO blanks produces dusts which can irritate the eyes, skin, and airways. As such, it is essential to ensure the protective gear at your workplace is in perfect working order.
- pro3dure medical GmbH is not liable for any damages caused by improper application of the material. To be used by trained specialist personnel for the purpose indicated only.

### 7. Label symbols

**Manufacturer:**  
pro3dure medical GmbH  
Am Burgberg 13 · 58642 Iserlohn, Germany  
Phone: +49 (0)2374 920050-0

**Distributor (US):**  
pro3dure medical LLC  
9825 Valley View Road · Eden Prairie, MN 55344  
Phone: +1 (952)-426-1928

info@pro3dure.com · www.pro3dure.com · Made in Germany