**1. Intended use / Indication**  
otosil® protect 50 | MSI from pro3dure is a 1:1 Silicone material for manufacturing of hearing protections.

**2. Contraindication**

otosil® protect 50 | MSI is contraindicated ...  
1. ... if a patient is known to be allergic to any of the ingredients.  
2. ... for every application that is not part of the indication (see above).

**3. Product Description**

The otosil® protect 50 | MSI includes a lactam-based, biomimetic solution for multi-species inhibition\*. This reduces the formation of undesirable biofilms on audio appliances without a microcidal effect.  
\*c. albicans / s. mutans

**4. Patient Target Group**

Persons for whom an earmold is to be created.

**5. Intended Users**

Earmold laboratory, ENT specialist, hearing aid acoustician

**6. Requirements**

Use by specialists only.

**7. Material**

Contains: Functionalized polydimethylsiloxanes, crosslinkers, modifiers and lactam.

**8. Manufacturing process (fig. 1-6)**

**Mixing and dosing**

- Place the cartridge into the dispensing gun.
- Unscrew cap. In order to ensure a proper flow and mixing ratio from both orifices extrude carefully a small amount of the material.
- Insert the static mixer into the guided grooves of cartridge. Rotate the mixing cannula anti-clockwise to a stop. If necessary, mount an Intra-tip. For reducing forces during dispensing, static mixers of a larger diameter (9 mm) can be used alternatively. Dispenser is now ready to mix the 2 component earmold silicone in any amount needed.

**Application (PnP-method)**

- With the otosil® protect 50 | MSI material all commercially available materials for negative forms like plasters, gels or 3D-printing resins can be used for the manufacturing of the negative form.
- If necessary, coat the negative form with a separation liquid (e.g. alginate based for plaster). 3D-printed cast forms have to be cleaned carefully before use in order to avoid inhibition of the silicone reaction. Inject the otosil® protect 50 | MSI material slowly and bubblefree into the generated negative form. In case of narrow structures (e.g. auditory canals), we recommend using a thin mixing canula. Put the negative mold with the injected material into a pressure pot for about 60 min at 122 °F (50 °C).

**Final surface treatment and lacquering**

- The silicone earmold can now be shaped by using special cutters and grinding sleeves. Roughen the surface with corundum paper (grit 180) before varnishing. Clean the surface and coat it with a silicone lacquer according to the special instructions of the lacquer.


**9. Important working hints**

In case of the generative fabricated castforms a careful preparation as well as a thorough cleaning will ensure better and more reliable results. Ensure that all traces of residual resin is removed from the negative form. The cured earmold material otosil® protect 50 | MSI is chemically inert. Avoid spots on clothing. We recommend standard gloves made of nitrile or polyethylene.

**10. Notice**

Undesired reactions of the immune system like allergies, irritations can't be absolutely excluded. In case of doubt, we recommend to contact your consulting physician and make an allergy test before the application of the material. According to the EU Medical Devices Regulation, users/patients are obliged to report serious events with a medical device to the manufacturer and to the competent authority of the country in which they occurred. Dispose of contents/container according to official regulations.

**Earmold silicone  
protect 50 | MSI 50 Shore**  
Instructions for use

**Physical properties\*:**

**otosil® protect 50 | MSI**

- Mixing volume: 50 ml
- Mixing ratio: 1:1
- Setting time (50 °C): ≤ 15 min. \*\*
- Processing time (23 °C): ≥ 4 min. \*\*
- Tensile strength / MPa (ISO 37:2017): ≥ 4.0 \*\*
- Elongation at break / % (ISO 37:2017): ≥ 200 \*\*
- Tear resistance / N/mm (DIN ISO 34:2016): ≥ 12 \*\*
- Shore hardness A (DIN ISO 48-4:2021): 50±5 \*\*

**Ordering information:**

**otosil® protect 50 | MSI**

**Standard package:**

Cartridge  
8 x 48 ml  
(à 2 x 24 ml A + B)

**neon-yellow**  
REF: A0034300

**neon-orange**  
REF: A0034301

**neon-pink**  
REF: A0034302

**Cannula:**  
A System  
5,4 mm  
transparent  
50 pieces

REF: A4020142

**Dispenser:**  
DMA50  
REF: A4020128

\* These data come from measurements of a representative sample that was determined as part of our quality assurance.  
\*\* According to internal design and requirements specifications.

This side up.

Temperature limit

CE mark

Do not use if package is damaged.

Consult instructions for use.

Keep dry.

Rx only

QTY: 1EA

Medical Device

Use by date

Catalogue number

Lot number

Date of manufacture

**Manufacturer:**  
pro3dure medical GmbH  
Am Burgberg 13 · 58642 Iserlohn, Germany  
Phone: +49 (0)2374 920050-0

**Distributor (US):**  
pro3dure medical LLC  
9825 Valley View Road · Eden Prairie, MN 55344  
Phone: 952-426-1928

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