PRO3dent[®] splint | MSI

Instructions for use
Gebrauchsanweisung
Mode d'emploi
Instrucciones de uso
Istruzioni per l'uso

pro3dure
medical

Materials are our DNA!

Milling disc

1. Intended use

PRO3dent[®] splint I MSI milling discs for the fabrication of functional therapeutic occlusal splints.

2. Indication

The **PRO3dent[®] splint I MSI** blank is a disc intended for the fabrication of dental appliances such as mouthguards, nightguards, bruxism and TMJ splint appliances, and bite splints using CAD/CAM milling machines.

3. Contraindication

PRO3dent[®] splint I MSI is contraindicated ...

- ... if a patient is known to be allergic to any of the ingredients.
- ... for bruxism splints and splints with wall thickness ≤ 0.9 mm.
- ... for denture bases.
- ... for every application that is not part of the indication (see above).

4. Product Description

The **PRO3dent[®] splint I MSI** includes a lactam-based, biomimetic solution for multi-species inhibition*. This reduces the formation of undesirable biofilms on orthodontic appliances without a microcidal effect. *c. albicans / s. mutans

5. Patient Target Group

Persons being treated in the context of a dental procedure.

6. Intended Users

Dentist, dental technicians

7. Processing

- PRO3dent[®] splint I MSI** milling blanks should be used only for the purpose of CAD/CAM production of dental occlusal splints. They can be used in all common CAD/CAM systems (fig. 1).
- Suitable cross-cut carbide cutters are to be used exclusively for the processing of the **PRO3dent[®] splint I MSI** material.
- PRO3dent[®] splint I MSI** occlusal splints are designed and milled by CAD/CAM technology and qualified staff. After finishing, the thickness in the occlusal region of the splint must not be less than 0,9 mm. For aesthetic reasons, a labial reduction to 0,8 mm is permissible. Any corners and edges should be rounded off.
- PRO3dent[®] splint I MSI** splints can be removed from the milling blank by using suitable cross-cut carbide cutters or cutting discs suitable for composites. Carefully cut through the retaining strips without pressure.
- In order to avoid plaque accumulation subsequent polishing of the splints is essential. Please pre-polish with suitable silicone

polishers and goat's hair brushes. The high-luster polishing should be done with a corresponding polishing compound (fig. 2). As far as possible, please avoid heat generation during the polishing and finishing of the **PRO3dent[®] splint I MSI** splint. This guarantees the optimal fitting of the splint.

- Please note:** Upon completion, use an ultrasonic bath with water at a maximum temperature of approx. 35 °C for cleaning. Cleaning solutions are not recommended.

8. Milling parameters

Roughing:

- 2 mm ball radius cutter – single cutter
- Speed: 22,000 rpm
- Radial path distance: 0.6 mm
- Z depth increase: 0.6 mm
- Feed: 1,600 mm/min

Finishing:

- 2 mm ball radius cutter – single cutter
- Speed: 22,000 rpm
- Radial path distance: 0.1 mm
- Z depth increase: 0.1 mm
- Feed: 1,600 mm/min

Rest roughing:


- 1 mm ball radius cutter – single cutter
- Speed: 28,000 rpm
- Radial path distance: 0.1 mm
- Z depth increase: 0.1 mm
- Feed: 1,200 mm/min

The irrigation nozzle on the spindle should be aligned with the tool tip. In this occasion it should be noted that a uniform programming length is respected or the length of the 2 mm tool is estimated. The cooling medium air is sufficient.

9. Notes

- Incorrect use of milling tools, polishing brushes, steam cleaner and water baths can overheat the material exposed and thus adversely affect the properties. Deviations from the described manufacturing process can lead to different mechanical and optical properties of the material.
- The lot number and the best before date are indicated on each **PRO3dent[®] splint I MSI** packaging. In case of claims please always indicate the lot number of the product.

UMDNS 16-697

- Do not use the product after expiration date.
- Undesirable biological reactions (e.g. allergies to material components) may occur in individual cases. Clarify the use of **PRO3dent[®] splint I MSI** splints with your dentist or orthodontist in case of a known allergic response to components of **PRO3dent[®] splint I MSI**.
- Clean **PRO3dent[®] splint I MSI** splints under cold running water with a toothbrush. Disinfect the splint with non-alcoholic detergents before first use.
- According to the EU Medical Devices Regulation, users/patients are obliged to report serious events with a medical device to the manufacturer and to the competent authority of the country in which they occurred.

10. Safety advice

- Record the BATCH no. for every procedure requiring identification of the material.
- Storage:** No special measures are required.
- Disposal:** Dental objects made of **PRO3dent[®] splint I MSI** are not water soluble, are inert, and do not present a hazard for ground water. As such, they can be disposed of as domestic waste without any special precautions if not stated otherwise in your local disposal guidelines.
- Warranty:** Our recommendations concerning the application technique, regardless of whether they are communicated in writing, orally, or by means of practical instructions, are based on our own experiences and tests. As such, they are intended as guidelines only. We are continually striving to improve our products. Consequently, we reserve the right to make changes to their design and composition.
- Hazard Information:** Processing of **PRO3dent[®] splint I MSI** blanks produces dusts which can irritate the eyes, skin, and airways. As such, it is essential to ensure the protective gear at your workplace is in perfect working order.
- pro3dure** is not liable for any damages caused by improper application of the material. To be used by trained specialist personnel for the purpose indicated only.

Physical properties */
Physikalische Eigenschaften */
Propriétés physiques */
Propiedades físicas */
Proprietà fisiche *:

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
Distributor (US):
pro3dure medical LLC
9825 Valley View Road · Eden Prairie, MN 55344
Phone: +1 (952)-426-1928

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