pro 3d ure
medical

D E N T A L

25 | 26

HarzRock

p3d

The Harz of Rock 'n' Roll

pro3dure
medical

DENTAL

25 | 26


Inhalt | Contents

4	Einleitung Introduction
6	printodent® NEWS
8	MSI® Technologie Technology
10	printodent® GR-10 guide _ GR-10 guide MSI _ GR-10.1 guide MSI
16	printodent® GR-11.1 tray
18	printodent® GR-12 cast
19	printodent® GR-13 model _ GR-13.1 model _ GR-13.2 model _ GR-13.3 model release _ GR-13.4 model NEW
24	printodent® GR-14.2 denture HI _ GR-14.2 denture HI MSI _ GR-14.3 denture HI MSI NEW
29	printodent® GR-15 gingiva
30	printodent® GR-16 Xray
31	printodent® GR-17 temporary _ GR-17.1 temporary It _ GR-17.2 temporary It MSI 
36	printodent® GR-18.3 IB NEW
38	printodent® GR-19.1 OA _ GR-19.1 OA MSI
40	printodent® GR-20 MJF
42	printodent® GR-21.1 Try-In
44	printodent® GR-22 flex
50	printodent® GR-23 mouthguard NEW
52	PRO3dent® Familie Family 
56	THERMEO® Familie Family
64	ISM-1 MSI mouthguard 
65	printisol Isoliermittel printisol Insulating agent NEW
66	Zubehör Accessories
68	Geräte (3D-Drucker, Reinigung, Nachhärtung) Devices (3D printer, cleaning, post curing)
78	Validierung Validation

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Future innovations are marked with the **innovation icon**.

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Materials are our **DNA!** ... is the guiding principle of pro3dure medical. Just like in our name, the 3d is the central pivot of all our activities and has been since 2013.

The pro3dure team's outstanding passion for materials used in digital processes is reflected in a special way in this year's IDS motto "The Harz of Rock 'n' Roll". It is therefore a matter close to our hearts to passionately develop and advance digital processes and the associated materials. That is why we are already showing you our future innovations, which you can identify by the **innovation icon** in the brochure. We are delighted that the approval has finally been given for the construction of our new competence center for medical 3D printing. One of the world's most modern production facilities for 3D printing materials, a research and development landscape covering all areas of chemistry and a training center for digital processes for the manufacture of medical products are now being built on almost 5000 m² in the Dortmund Technology Park. You can follow all of this in our construction diary: www.pro3dure.com/en/construction-diary. We have set ourselves the goal of learning from and with you in order to innovatively advance the manufacture of dental applications in all areas. Trusting cooperation and proximity to our customers and partners is of the highest of importance to us in our daily work. Let's take the next steps together!

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printo[®] dent

MSI[®] technology

printo[®] GR-10 guide
 printo[®] GR-10 guide | MSI
 printo[®] GR-10.1 guide | MSI
 printo[®] GR-11.1 tray
 printo[®] GR-12 cast
 printo[®] GR-13 model
 printo[®] GR-13.1 model
 printo[®] GR-13.2 model
 printo[®] GR-13.3 model release


printo[®] GR-13.4 model **NEW**
 printo[®] GR-14.2 denture HI
 printo[®] GR-14.2 denture HI | MSI
 printo[®] GR-14.3 denture HI | MSI **NEW**
 printo[®] GR-15 gingiva
 printo[®] GR-16 Xray
 printo[®] GR-17 temporary
 printo[®] GR-17.1 temporary It
 printo[®] GR-17.2 temporary It | MSI

printo[®] GR-18.3 IB **NEW**
 printo[®] GR-19.1 OA
 printo[®] GR-19.1 OA | MSI
 printo[®] GR-20 MJF
 printo[®] GR-21.1 Try-In
 printo[®] GR-22 flex
 printo[®] GR-23 mouthguard **NEW**

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
The world's largest portfolio of 3D printing resins for dental applications becomes even more comprehensive. **printo[®]** resins with the worldwide unique **MSI[®] technology** are especially marked by the **MSI[®] Logo**.

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
Communication is everything and represents one of the most fundamental mechanisms of intra- and interspecies interactions. Communication can take place in many different ways. One form of human communication is the acoustic transmission of signals to transfer information. In addition to acoustic communication, visual communication also plays a major role in the animal kingdom. In an impressive way, this can be observed in marine-living cuttlefish (octopus). These visually communicating organisms are not only able to change their skin color and surface so that they are perfectly camouflaged. They are also able to move platelets in the chromatophores (specialized skin cells) to reflect specific patterns of polarized radiation. It is now thought that they have a language via these different skin changes, consisting of over 50 vocabularies, and can even form sentences in their specific sequence. Both acoustic and visual communication require that the sender and receiver have highly complex and specialized organs such as ears and eyes.

But what happens when, as in the case of bacteria, the organism cannot possess these organs because it consists of only a single cell? Although bacteria are only single-celled, they can communicate with conspecifics and bacteria of other species. This mechanism was first described in symbiotic bioluminescent bacteria by Nealson et. al. in 1970 and termed **Quorum Sensing** by Fuqua in 1994 and represents a biochemical language. The vocabulary of this language are so-called autoinducers, which are formed by the bacteria themselves and released into the space surrounding. The more bacteria exist in a spatial proximity to each other, the more of these signal molecules are present in this environment. When the concentration of the messenger molecules reaches a certain threshold value (the quorum), the autoinducers flow back into the bacteria, switch on some genes, switch off others and thereby abruptly change the behavior of the organisms. Sometimes the microorganisms form a biofilm that protects them from antibiotics and disinfectants; mouth germs then form harmful dental plaque. In other cases, they produce toxins and decomposing enzymes.

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(<https://onlinelibrary.wiley.com/doi/10.1111/apm.13510>)

Prof. Gordon Ramage | Glasgow Caledonian University:

"The work carried out at Glasgow Caledonian University has shown that **MSI[®] technology** is able to inhibit and prevent biofilm formation when incorporated into dental materials. This is most obvious when we investigate yeasts." The Ramage lab has recently reported in an international journal that the parent technology for **MSI[®]** is highly biofilm active.

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Bacteria have a fascinating and diverse social life. They exhibit coordinated group behavior controlled by **quorum-sensing systems** that detect the density of other bacteria in their environment. The term quorum dates back to the Roman Empire and referred to the lowest number of members required for a vote in the Senate. Similar to the Roman Senate, bacteria vote their behavior as soon as a certain number of other bacteria are in their vicinity. A key example of such group behavior is biofilm formation (e. g. plaque), in which communities of cells attach to a surface and form a polymer shell that protects the bacteria from external influences.

multi-species inhibition

MSI[®] technology Inspired by Nature | Addicted to Function

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Prof. Richard Hammond | Penrhos Bio Chief Executive Officer:
“We are thrilled to be working with our colleagues **pro3dure medical** to combine the power of **Remora[™]** technology with their world leading expertise in material development. The result is a completely novel way of inhibiting the formation of harmful multispecies biofilms in their products. **Remora[™]** is inspired by the natural biology of certain seaweeds that are able to block communication between microorganisms. The result is a surface that remains clean and free of biofilms. Backed by Unilever, we have proven efficacy in a wide range of applications and we are really excited to bring this benefit into dental markets through the innovation of **pro3dure**, an ideal partnership.“

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Penrhos Bio, Unilever, **pro3dure medical GmbH** – Winner of Scotland's Life Sciences Annual Awards 2024 Glasgow (Category Innovative Collaboration). A great evening with our friends and partners was highlighted with a big win for our collaboration! In addition, Glasgow Caledonian University, **pro3dure** and Penrhos Bio were awarded the National Biofilms Innovation Centre (NBIC) grant.

The formation of biofilms is therefore a significant medical problem (e. g. caries formation). If the communication of microbes could now be interfered with by means of suitable molecules, bacterial infections could also be contained. Unlike antibiotics, the germs would not be killed, but cut off from their information channels. For this reason, the organisms cannot develop resistance to the interfering molecules, as is known from antibiotic resistance. (Fig. 4) In nature, biofilms are not observed on all moist surfaces. There are species of algae that are virtually free of biofilms. For example, the red alga *Delisea* produces a variety of natural compounds (furanones) that appear to block quorum sensing by bacteria. Inspired by nature, **Unilever / Penrhos Bio** has developed a biomimetic technology based on these natural compounds to combat biofilms (called **Remora[™]**). Research results also document that the bacteriostatic effect works not only on bacteria but also on yeasts, fungi and algae. Thus, a technology is available that can effectively combat „multi-species“ biofilms.

In close cooperation between **Penrhos Bio** and **pro3dure medical**, this multi-species inhibition (**MSI[®] technology**) was realized for the first time using natural substances in 3D printing materials and other dental materials. In doing so, **pro3dure medical** is taking an innovative step towards functional materials in the further development of its dental materials in accordance with the **MSI[®] guiding principle** "Inspired by nature, addicted to function". Come and find out more details about **MSI[®] technology** and see numerous new, functional dental materials ranging from milling blanks and denture relining materials to 3D printing resins.

Abb. 3: printo^{dent}[®] Harze mit MSI[®] Technologie
Fig. 3: printo^{dent}[®] resins with MSI[®] technology

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pro3dure's 3D printing resins of the **printo[®]dent[®] GR-10 series** are biocompatible resins for the fabrication of splints and implant drill guides. In addition to its exceptional elongation for this class of material, **printo[®]dent[®] GR-10 guide** is also characterized by its high transparency and UV stability. This makes it easy to insert sleeves into the guide and gives you optimum control during the operation. The material is now also available in the **printo[®]dent[®] GR-10 guide | MSI** version as a functional 3D printing resin with **MSI[®] technology**, which significantly reduces biofilm and plaque build-up on surfaces. **MSI[®]** is based on the patented mode of action of special natural compounds that disrupt the communication of bacteria (quorum sensing). Inspired by these natural compounds, which are produced by the red alga *Delisea*, a sustainable approach to preventing biofilms has been found that, unlike antibiotics, for example, does not cause bacterial resistance.

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MSI[®] technology in combination with a new polymer matrix has been realized in the **printo[®]dent[®] GR-10.1 guide | MSI** material. In addition to the advantages of reduced biofilm adhesion, the new formulation results in outstanding mechanical properties such as fracture toughness and fracture work. Implant drill guides made of **printo[®]dent[®] GR-10.1 guide | MSI** meet the highest stability requirements and have a fracture toughness and fracture work that are many times higher than those of their market competitors (Fig. 8/9, page 14). This minimizes the risk of splint fracture and increases clinical safety. Due to its novel formulation, it has been possible to dispense with teratogenic components, thus significantly increasing body compatibility. All materials of the **printo[®]dent[®] GR-10 series** can be sterilized in autoclaves using validated processes (Fig. 6).

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pro3dure's printo[®] dent GR-11 tray materials are biocompatible 3D-printing resins for the additive manufacturing of customized impression and functional impression trays. Due to their high reactivity, buildstyles with up to 300 µm Z-height can be achieved to be efficient in the lab. With the new development of pro3dure's printo[®] dent GR-11.1 tray, the biocompatibility and mechanical stability of the built objects have been taken to another level. The new formulation does not have any reproduction-toxic labelling, thus documenting the special biocompatibility for users and patients. Furthermore, pro3dure's printo[®] dent GR-11.1 tray has an impact strength that is approximately 50 % higher than any competitor product (Fig. 10). Accordingly, impression trays can be removed from the patient's mouth with higher safety.

pro3dure's printo[®] dent GR-11 tray materials are biocompatible 3D printing resins for the additive manufacturing of customized impression and functional impression trays. Due to their high reactivity, buildstyles with up to 300 µm Z-height can be achieved to be efficient in the lab. With the new development of pro3dure's printo[®] dent GR-11.1 tray, the biocompatibility and mechanical stability of the built objects have been taken to another level. The new formulation does not have any reproduction-toxic labelling, thus documenting the special biocompatibility for users and patients. Furthermore, pro3dure's printo[®] dent GR-11.1 tray has an impact strength that is approximately 50 % higher than any competitor product (Fig. 10). Accordingly, impression trays can be removed from the patient's mouth with higher safety.

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1. Preheat oven to 600 °C.
2. Set the muffle.
3. Holding time¹ 45 min.

¹ Holding time/casting process according to manufacturer's instructions

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- rückstandsfrei ausbrennbar
- hohe Präzision, feinste Details
- hohe Stabilität
- schnelle, effiziente Reinigung
- reparierbar mit LCR-1
- residue-free burn-out
- high precision, finest details
- high stability
- fast, efficient cleaning
- repairable with LCR-1

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pro3dure's printo[®]dent GR-12 cast material is a residue-free burn-out 3D printing resin that can be used in combination with a variety of investment materials for a wide range of applications. Thanks to the optimized setting of the curing depth, even the finest details can be printed precisely with the material. The high stability of the material in the green compact mode even allows a direct FIT check. Components made of **printo[®]dent GR-12 cast** can be cleaned quickly and efficiently thanks to its low viscosity. If, despite their high stability, **printo[®]dent GR-12 cast** objects need to be repaired or augmented, we recommend the associated UV-curable repair gel **LCR-1** from **pro3dure** (see page 68).

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“3D printing with **printo[®]dent** in your own medical office is the perfect complement to our intraoral scanner.”

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Fig. 14: Model of the same situation printed with **printo[®]dent GR-13 model**

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pro3dure's printo[®]dent GR-13 model is a 3D printing resin for the generative fabrication of dental models. The material is available as **printo[®]dent GR-13 model** for restorative dental models and as **printo[®]dent GR-13.1** / **printo[®]dent GR-13.2 model** for orthodontic applications. All resins are characterized by high dimensional stability and the optimized look and feel for dental applications. Due to colorant and opacity the important features like margin lines are becoming perfectly visual. Moreover, the implementation of 3D printed models into your digital workflow is supported. Dental models from the **GR-13 series** can be used effectively in thermoforming processes due to their thermal stability.

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Fig. 15: Release forces of different model materials in conjunction with an orthodontic powder/liquid system (ISO 22112), N ***

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printo[®]dent GR-13.3 model release was specially developed for orthodontic applications and eliminates the need for additional release agents when using 2-component powder/liquid resins with 3D-printed dental models. Through pro3dure's new "release" technology, it is possible for the first time to provide a material that has an integrated separation effect. This is not due to the addition of fillers such as waxes, which can segregate from the resin matrix and lead to unstable builds or defective components. Furthermore, due to its low viscosity, components made of printo[®]dent GR-13.3 model release are easy to clean, thus enabling an efficient finishing process. It is suitable for many orthodontic powder/liquid systems[°].

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The new **printo[®]dent GR-13.4 model stone** is the first model material in its class to give you the real “plaster feeling”. Thanks to its plaster-like feel in combination with high density, surface hardness and strength, you can enjoy the advantages of the “analog” and “digital” model world at the same time. Dental models made of **printo[®]dent GR-13.4 model stone** have a matte surface that offers optimum contrast even with the finest details. Furthermore, **printo[®]dent GR-13.4 model stone** is compatible with both **printisol** and conventional, water-based insulating materials. Another advantage of **printo[®]dent GR-13.4 model stone** is its high temperature stability, which is also beneficial when working with Bunsen burners and soldering irons.

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printo[®]dent GR-13.4 model stone
1 kg Flasche / bottle

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
These data result from measurements of a representative sample, which were determined within the scope of our quality assurance.

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printo[®]dent GR-13 Serie / series

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
a) Bauplattform mit Arbeitsmodell (inklusive separater Gingiva-Masken) aus **printo[®]dent GR-13.3 model release**. / Building platform with working model (including separate gingival masks) made of **printo[®]dent GR-13.3 model release**.

b) Gedrucktes Dental-Modell aus **printo[®]dent GR-13 model** (beige) und Arbeitsmodell aus **printo[®]dent GR-13.3 model release** (weiß-opak) in Kombination mit 3D-gedruckter bzw. gefräster Basis und gefrästen Prothesenzähnen. / Printed dental model from **printo[®]dent GR-13 model** (beige) and working model from **printo[®]dent GR-13.3 model release** (white-opaque) in combination with 3D-printed or milled base and milled denture teeth.

c) 3D-gedruckte Modelle aus **printo[®]dent GR-13 model** (beige) als Design-Mastermodell. **printo[®]dent GR-13.3 model release** (weiß-opak) als Modell zum Befestigen der Locatoren. / 3D-printed models made of **printo[®]dent GR-13 model** (beige) as design master model. **printo[®]dent GR-13.3 model release** (white-opaque) as model for attaching the locators.

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Fig. 16 (a-e): Use of different model materials (**printo[®]dent GR-13 model** and **13.3 model release**) in the digital workflow.

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Combination of different
model materials
in the digital workflow

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
Materials are our DNA! 23

printo[®]dent GR-14.2 denture HI

high accuracy of fit
hohe Passgenauigkeit

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- hohe Passgenauigkeit
- „high-impact“ Material
- validierte Verfahren
- high accuracy of fit
- „high-impact“ material
- validated processes

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pro3dure offers some of the finest denture base and tooth materials on the market today. By seamlessly blending strength and aesthetics, they have established their products as a top choice within the dental industry. The high-impact material ensures patients receive durable, long-lasting dentures. In my dental laboratory this base is the number one base used on every denture.

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The **printo[®]dent GR-14 denture** materials are 3D printing resins for the generative production of denture bases. The **printo[®]dent GR-14.2 denture HI** meets all dental requirements with e.g. a water absorption (< 25 µg/mm³) and water solubility (< 1.5 µg/mm³) according to ISO 20795-1:2013. Furthermore, the low shrinkage (approx. < 1 %) compared to traditionally used PMMA-based materials guarantees a high accuracy of fit. **printo[®]dent GR-14.2 denture HI** is a 3D printing resin with „high-impact“ properties characterized by low viscosity. This allows construction processes to be carried out efficiently and components to be cleaned easily. As a world first, **printo[®]dent GR-14.2 denture HI | MSI**, a biomimetic 3D printing resin for the fabrication of dentures, was presented for the first time at IDS 2023. This resin offers all the advantages of **printo[®]dent GR-14.2 denture HI**, which has been tried and tested thousands of times, and also reduces biofilms from adhering to the surface of the denture in a natural way. This functional material thus copies a natural process, such as that used by the red alga *Delisea*. The communication of bacteria, the so-called „quorum sensing“, is disturbed by a component anchored in **printo[®]dent GR-14.2 denture HI | MSI**, so that the adhesion of bacteria is significantly reduced. Moreover, this effect is not due to antibacterial additives, which may lead to bacterial resistance. As a result, improved oral hygiene and less irritation of the soft tissue are achieved.

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- exceptional mechanical properties
- **MSI[®] technology** for the reduction of biofilms
- high accuracy of fit
- better color stability
- high stability after storage in water
- "high-impact" material
- TPO-free

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The result of our continuous development activities is the new **printo[®]dent GR-14.3 denture HI | MSI**, which is characterized by exceptional mechanical properties. Compared to its market competitors, the mechanical properties remain significantly more stable even after storage in water (Fig. 22). As a further member of the **printo[®]dent** family, the **printo[®]dent GR-14.3 denture HI | MSI** is also equipped with **MSI[®] technology** to reduce biofilm formation. In combination with outstanding chemical-physical properties, dentures made of this material have a thermo effect optimized for body temperature, which generates an extremely high level of wearer comfort. All materials of the **printo[®]dent GR-14 denture series** can of course be combined with **printo[®]dent GR-17.1 temporary It** to produce cost-effective and highly aesthetic full dentures. In addition, validated procedures are available for both materials in combination with denture teeth from the **VITA VIONIC VIGO**, **CANDULOR TCR[®] / NFC[®]** and **Ruthinium ACRYSMART** systems. Complete integration in **exocad**'s CAD software rounds off the digital workflow.

* CANDULOR Bonartic[®] II NFC[®], Bonartic[®] TCR, BonSelect TCR, Condyloform[®] II NFC[®], PhysioSelect[®] TCR, PhysioSet[®] TCR

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"Biofilm inhibition with the Remora™ brand technology results in cleaner and healthier dental products and appliances which assist in maintaining oral health and stability."

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Fig. 22: Maximum stress intensity factor – mechanical change after water storage, % *

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pro3dure's printodent® GR-15 gingiva is a 3D printing resin for the production of flexible gingiva masks. Its natural-looking gingival color combined with a gingival-like flexibility allow optimal dental work with a view to "red and white aesthetics". The printodent® GR-15 gingiva is the perfect complement to pro3dure's printodent® GR-13 model to fabricate implant models based on biometric data.

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"To ensure a harmonious transition from the crown to the gum, printodent® GR-15 gingiva is excellent."

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- natural coloration
- gingiva-like flexibility
- perfect complement to printodent® GR-13 model

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"We print a radiopaque scanning appliance from **printodent[®] GR-16 Xray** with each of our end of the line digital dentures and All-On-X cases. This allows us to distinctly differentiate our lab with our digital workflow and deliver a better prostheses to our dentists."

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pro3dure's printo[®]dent GR-16 Xray is the first radiopaque 3D printing resin for the generative fabrication of scanning templates. On the basis of a stabilization technology developed by **pro3dure** for filled materials and a suitable selection of radiopaque additives, the construction of scanning templates can now be integrated into your "digital workflow" efficiently.

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Fig. 24:
Radiopacity of **printo[®]dent GR-16 Xray** in comparison to aluminum stairs (ISO 4049).

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- biocompatible materials
- very good mechanical properties
- double indication
- 7 colors - 2 sizes
- low viscosity
- easy to clean
- radiopaque

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pro3dure's GR-17 series resins are biocompatible materials for the additive manufacturing of temporary crowns and bridges. The **printodent[®] GR-17 temporary** is indicated for temporary restorations in the anterior region. **printodent[®] GR-17.1 temporary It** has exceptional mechanical properties compared to conventional temporary materials with a Flexural modulus > 5000 MPa (according to ISO 178**) and a Flexural strength of 112 MPa (ISO 4049). Moreover **printodent[®] GR-17.1 temporary It** exhibits excellent tissue compatibility and can therefore be used ideally for longterm temporaries¹. As the first material in its class worldwide, **printodent[®] GR-17.1 temporary It** received FDA approval for two indications: temporaries and denture teeth. This expands the possibilities in your laboratory and makes the use of this material extremely efficient and economical. Both materials are available in tooth shades A1, A2, A3, A3.5, B1, B2 and bleach.

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
Fig. 26: Radiopacity of printo[®]
GR-17.1 temporary It in comparison
to aluminum stairs (ISO 4049).

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For the temporary restoration after
implantation, I found the optimal
solution for my digital workflow
with printo[®] GR-17.1 temporary It.

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Scientific studies¹ prove the compatibility with bredent visio.link as well as the
optimal adhesive bond with printo[®] GR-17.1 temporary It.

¹Lankes V, Reymus M, Liebermann A, Stawarczyk B. Bond
strength between temporary 3D printable resin and con-
ventional resin composite: influence of cleaning methods
and air-abrasion parameters. Clin Oral Invest. 2023
Jan;27(1):31-43. doi: 10.1007/s00784-022-04800-7

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- highest mechanical properties in its class
- TPO-free
- filler content 60 m %
- biocompatible
- double indication
- 3 colors – 2 sizes
- easy to clean
- radiopaque
- **MSI[®] technology** for the reduction of biofilms

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
pro3dure's **GR-17 series** resins are biocompatible materials for the additive manufacturing of temporary crowns and bridges. Thanks to its high filler content of 60 m %, **printo[®]dent GR-17.2 temporary It | MSI** has the highest mechanical properties in its class and is easy to clean. **printo[®]dent GR-17.2 temporary It | MSI** has outstanding mechanical properties and excellent soft tissue compatibility compared to conventional temporary materials with a flexural modulus > 7000 MPa (according to ISO 178 **) and a flexural strength of 123 MPa (ISO 4049). The additional implementation of **MSI[®] technology** and the associated reduction of biofilms make the material ideal for long-term temporary restorations in combination with implants.

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
pro3dure's printo[®]dent[®] GR-18.3 IB is a biocompatible resin for the additive manufacturing of orthodontic indirect bonding trays. It is the enhanced development of the tried-and-tested printo[®]dent[®] GR-18.2 IB in terms of biocompatibility and process reliability. The material is flexible, transparent and has even greater tear resistance. This allows the orthodontist to place all brackets in one controlled, quick and precise step. Biocompatibility was a main focus during further development. The new printo[®]dent[®] GR-18.3 IB is now TPO-free. Although the new formulation has further increased the speed and stability of the construction process, printo[®]dent[®] GR-18.3 IB is still easy to clean. The inhibition layer during post-curing with low-power post-curing systems is still minimal, thus avoiding sticky component surfaces and ensuring greater clinical safety.

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

pro3dure's printo[®] dent GR-19.1 OA | MSI is the first biomimetic 3D printing resin for the generative production of hard-flexible orthodontic splints. Thanks to the unique **MSI[®] technology**, splints made of **printo[®] dent GR-19.1 OA | MSI** exhibit reduced biofilm adhesion. The effect, which is copied from nature, is based on components that sustainably inhibit the communication of bacteria. Unlike antibiotics, for example, no resistance is formed in this way. The material also has a body temperature-induced thermo effect and thus ensures maximum wearing comfort in the mouth.

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"I use **printo[®] dent GR-19.1 | MSI** because I make good use of milling capacity in my laboratory. I have been offering the 3D printed **MSI[®]** splints for 1 year and have an almost 100% success rate. The **MSI[®] technology** for reducing biofilms has been very well received by my customers."

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pro3dure's printo[®]dent[®] GR-20 MJF is a biocompatible and hemocompatible tested (according to ISO 10993-4) 3D printing resin for special objects such as guides and splints.

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printo[®]dent GR-21.1 Try-In
1 kg Flasche / bottle

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printo[®]dent GR-21.1 Try-In is a biocompatible class 1 medical device and is suitable as base material for functional try-in of digitally fabricated denture bases. **printo[®]dent GR-21.1 Try-In** can be used to generatively fabricate dental moldings such as denture bases with individual tooth set-up in order to perform a functional try-in to check bite registration and occlusion. The material is available in tooth shades A1, A2 and A3. The newly developed **printo[®]dent GR-21.1 Try-In** is characterized by significantly improved chemical and physical values compared to its market competitors. For example, the impact strength (according to ISO 180/U) has been increased by 75 % compared to its predecessor. The use of a novel initiator system has further increased biocompatibility and thus clinical safety.

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
Fig. 29 (a-d): Try-in on an implant bridge with **printo[®]dent GR-21.1 Try-In**

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"It has long been my wish to be able to produce flexible splints without impressions. With **printodent[®] GR-22 flex**, we can now print our flexible splints in high quality. A significant step forward on my digital path."

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printodent[®] GR-22 flex – the innovation for 3D printing.

The number of treatments of occlusal disorders by means of splint therapy is constantly increasing and thus the importance of the "dental splint" for a wide variety of indications such as bruxism, bite elevations or mouth protection. The demands placed on splint materials by patients, practitioners and dental technicians are manifold and range from durability, aesthetics, wearing comfort and biocompatibility to processability and usability in digital processes. Based on biometric data, 3D printing processes can nowadays be used to produce complexly designed therapeutic splints efficiently. The latest member of the printodent[®] family was developed with these aspects in mind – the innovative, flexible splint material printodent[®] GR-22 flex with thermomemory effect, optimized for additive manufacturing. After all, who doesn't know it, the challenge of a "crooked impression" during conventional impression taking or imperfections during intraoral scanning that lead to "data holes" on the surface of a generated data set? A conventional, hard splint then wobbles on the plaster model and in the patient's mouth. The splint is subjected to greater stress at precisely these points. The risk of chipping or breaking is significantly increased. In addition, the wearing comfort for the patient is reduced. The innovative printodent[®] GR-22 flex compensates for such inaccuracies thanks to its elastic properties combined with a thermomemory effect. printodent[®] GR-22 flex becomes flexible at body temperature and thus adapts optimally to the tooth surface. This increases wearing comfort and minimizes the risk of fracture. printodent[®] GR-22 flex thus supports the satisfaction of dental technicians, practitioners and patients through higher safety, exceptional wearing comfort and fewer "recalls".

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be flexible

printodent® GR-22 flex is the flexible 3D printing resin for splints.

printodent® GR-22 flex gives you and your customers the flexibility you need! The development of the printodent® GR-22 flex was based on extensive clinical and material data from conventional, flexible splint materials. As a result, the printodent® GR-22 flex achieves increased flexibility in combination with a more pronounced thermomemory effect compared to its market competitors (Fig. 30), which is immediately noticeable when the splint is inserted in the mouth due to a higher pressure drop at the tooth. The wearing comfort for the patient is thus significantly increased. Would you also like to remain flexible in the selection and design of your "digital workflow"? A variety of validated processes based on a wide range of CAD softwares, 3D printers and post-treatment units are available for the printodent® GR-22 flex.

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printodent® GR-22 flex – a resin with outstanding properties.

printodent® GR 22 flex meets the highest medical technology standards. Dental splints made of printodent® GR 22 flex have outstanding mechanical properties compared to their market competitors. For example, the standard specifications for impact strength (ASTM_D256) and elongation at break (ISO 37) are not only achieved, but exceeded by more than 3.5 and 2.5 times respectively (Fig. 31 and Fig. 32). The new development printodent® GR 22 flex is based on an innovative polymer matrix which ensures that the solubility of polymer components (Fig. 33) is also significantly lower than that of the comparative products. This ensures outstanding biocompatibility in this class of materials.

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printodent® GR-22 flex is durable and biocompatible.

Characterized by its exceptional wearing comfort, printodent® GR-22 flex has been developed for durability and reliability. The stability of the material was proven in an extensive chewing simulation study conducted by LMU Munich¹. Although printodent® GR-22 flex has a much higher flexibility than its market competitor, no significant differences in terms of material loss (volumetric and vertical) could be found between them (Fig. 34 and Fig. 35). printodent® GR-22 flex thus increases patient comfort without reducing clinical benefit! Furthermore, tested reworking and cleaning recommendations are available to sustainably ensure the special properties of a printodent® GR-22 flex splint over its entire service life.

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Fig. 34: printodent® GR-22 flex (orange) and KeySplint Soft (blue) volumetric degradation after 120,000 cycles chewing simulation, mm³.¹

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Fig. 35: printodent® GR-22 flex (orange) and KeySplint Soft (blue) vertical degradation after 120,000 cycles chewing simulation, mm.¹

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printodent® GR-22 flex is the innovative resin for 3D printing of splints with unique impact strength combined with high elongation at break. Its special viscoplastic properties in combination with a thermomemory effect combine highest comfort with highest biocompatibility.

- high elongation at break
- unique impact strength
- highest wearing comfort (usability studies pro3dure)
- lowest solubility/high biocompatibility
- thermomemory effect
- highest impact strength in its class
- "snap on" splints available

printodent® GR-22 flex
1 kg Flasche / bottle
klar-transparent / clear-transparent
Art.-Nr. / item-no.: D1001701
hellblau / light-blue
Art.-Nr. / item-no.: D1001703

printodent® GR-22 flex
500 g Flasche / bottle
klar-transparent / clear-transparent
Art.-Nr. / item-no.: D1001702
hellblau / light-blue
Art.-Nr. / item-no.: D1001704

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
- high construction speed
- high process reliability
- high impact resistance
- high tear resistance
- optimal wearing comfort
- TPO-free
- available in 5 colors and clear-transparent

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pro3dure's printo[®]dent GR-23 mouthguard is a biocompatible 3D printing resin for the generative production of mouthguards. The material has been specially optimized for 3D printing and is characterized by high build speeds and process reliability. 3D-printed mouthguards made from printo[®]dent GR-23 mouthguard have a Shore D hardness of 40. Compared to other materials in this class, its mechanical values remain significantly more stable even after storage in water. The material is available in the colors white-opaque, blue-opaque, red-opaque, black-opaque, neon-orange-translucent and also in clear-transparent.

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
To complement our 3D printing resin portfolio **printodent[®]** and the high-tech **THERMEO[®]** resins, we bundle resins for a wide range of prosthetic applications under the umbrella brand **PRO3dent[®]**. For the digital and analog workflow, you will find technically outstanding products with the latest technology under the **PRO3dent[®]** brand. Optionally equipped with **MSI[®] technology**, **PRO3dent[®]** products offer a unique additional feature: multi-species inhibition significantly reduces biofilm accumulation and prevents plaque build-up.

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PRO3dent[®] splint | MSI is ideal for CAM-fabricated hard bite splints. Multi-species inhibition significantly reduces biofilm accumulation and prevents plaque build-up. Available in ocean-blue, in thicknesses of 16 mm and 20 mm (single packaging Ronde 98.5 mm with step).

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- hohe Biegefestigkeit
- keine Verfärbungen
- höchste Biokompatibilität
- Multi-Species-Inhibierung (**MSI[®]**)
- signifikante Reduktion von Biofilmen
- geringe Plaqueanlagerung
- high flexural strength
- no discoloration
- highest biocompatibility
- multi-species inhibition (**MSI[®]**)
- significant reduction of biofilms
- low plaque accumulation

PRO3dent [®] splint MSI	Ronde ozean-blau / blank ocean-blue	16 mm	Art.-Nr. / item-no.: D1025002
PRO3dent [®] splint MSI	Ronde ozean-blau / blank ocean-blue	20 mm	Art.-Nr. / item-no.: D1025003

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PRO3dent[®] base | MSI is used for the CAM production of individual denture bases in the digital workflow. Available in three colors and two thicknesses, the round blank with a diameter of 98.5 mm offers outstanding properties. The high flexural strength as well as the excellent bending modulus ensure optimum stability. Multi-species inhibition significantly reduces biofilm accumulation and prevents plaque build-up. Available in light pink, orange-pink and deep pink, in thicknesses of 25 mm and 30 mm (single packaging round 98.5 mm with step).

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- hohe Biegefestigkeit
 - keine Verfärbungen
 - höchste Biokompatibilität
 - Multi-Species-Inhibierung (MSI[®])
signifikante Reduktion von Biofilmen
 - geringe Plaqueanlagerung
- high flexural strength
 - no discoloration
 - highest biocompatibility
 - multi-species inhibition (MSI[®])
 - significant reduction of biofilms
 - low plaque accumulation

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The new **PRO3dent[®] flex base | MSI** discs offer multiple advantages, providing the option for a digital workflow in the creation of flexible denture bases. The high-performance acrylic is known for its exceptional wearing comfort, thanks to its thermomemory effect. Additionally, flexible denture bases are highly fracture-resistant and feature unique **MSI[®] technology**, which helps reduce biofilms and significantly enhances clinical safety. The material is available in three colors, carefully matched to the gums to seamlessly blend with the natural appearance of the mouth. With **PRO3dent[®] flex base | MSI**, metal retaining clasps, typically used in conventional dentures, are no longer necessary.

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- Multi-Species-Inhibierung (MSI[®])
signifikante Reduktion von Biofilmen
 - metallfrei und bruchstark
 - ästhetisch
 - höchster Tragekomfort durch Thermomemory-Effekt
- Multi-Species Inhibition (MSI[®])
significant reduction of biofilms
 - metal-free and unbreakable
 - aesthetic
 - highest wearing comfort due to thermomemory effect

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THERMEO® is flexible.

Who doesn't know it, "crooked impressions" during conventional impression taking or imperfections during intraoral scanning that lead to "data holes" on the surface of the generated data set? A conventional hard splint then wobbles on the plaster model and thus in the patient's mouth. The splint is subjected to greater stress at this point, the risk of fracture is significantly increased and the patient complains about a lack of wearing comfort. The unique THERMEO® splint material compensates for such inaccuracies with its patented thermomemory effect. The THERMEO® material becomes flexible at body temperature and adapts optimally to the tooth surface. This increases the wearing comfort of the splint and minimizes the risk of breakage. THERMEO® thus supports the satisfaction of the dental technician, dentist and patient through higher safety, exceptional wearing comfort and fewer "recalls".

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

"We are proud to be able to equip internationally successful world-class athletes in the field of sports dentistry with THERMEO® splints with thermomemory effect."

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THERMEO® is smart.

Smart materials/responsive materials are materials that react independently to changing environmental conditions such as temperature. This reaction is reversible. The high-tech plastics of the THERMEO® family belong to this group of smart materials and are characterized by a unique patented thermomemory effect. Compared to hard dental splints, inaccuracies during impression taking can thus be compensated with a THERMEO® splint. THERMEO® products can also be used to produce functional therapeutic bite splints. THERMEO® splints are characterized by high wearing comfort and are self-adjusting due to the thermomemory effect. THERMEO®s special properties combine the advantages of hard as well as soft splint materials.

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THERMEO® is durable.

The materials of the THERMEO® family, characterized by their exceptional wearing comfort, have been developed for durability. The stability of the material was documented in an extensive chewing simulation study conducted by LMU Munich*. Despite the higher flexibility and higher resilience of THERMEO® compared to its market competitors, comparable abrasion values were determined (Fig. 37 and 38). Thus, patient comfort is increased without reducing the clinical benefit. Furthermore, amine compounds and MMA were avoided in the development of THERMEO® in order to reduce the yellowing of the material and its potential for irritation and sensitization. Furthermore, tested cleaning recommendations for THERMEO® are available to the patient in order to sustainably ensure the special properties over the entire service life of the THERMEO® splint (see page 62).

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Fig. 37: THERMEO® (orange) and competitor A (red) vertical degradation after 120,000 cycles chewing simulation, mm.*

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Fig. 38: THERMEO® (orange) and competitor A (red) volumetric degradation after 120,000 cycles chewing simulation, mm³.*

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THERMEO® is comfortable.

THERMEO® meets the highest medical technology standards. For example, THERMEO® does not use phthalate-based plasticizers, but instead employs a patented flexibilization technology that ensures unique biocompatibility in this class of material. Compared to other suppliers, this results in a significantly greater pressure drop on the tooth after insertion of the splint (Fig. 40 and 41), which in turn increases wearing comfort for the patient. Furthermore, due to the technology used, the thermomemory effect and thus the resilience of the THERMEO® material is greater than that of competitors (Fig. 39). The THERMEO® splint returns to its original position faster and more precisely due to body temperature. THERMEO® can be processed conventionally using the casting or tamping-pressing technique or integrated into the digital workflow in the form of THERMEO® blanks. The blanks are available in thicknesses of 16, 20 and 25 mm. This allows tooth-supported and also soft tissue-supported splints to be manufactured efficiently. Both the powder/liquid system and the milling rounds are additionally available as THERMEO® SO in shades A1, A2 and A3 for the fabrication of "snap-on splints" (page 63).

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Fig. 39: THERMEO® (orange) and competitor A (red) elastic recovery after a 90° deformation depending on the temperature, 73.4 °F (23 °C) and 98.6 °F (37 °C).*

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Fig. 40: THERMEO® (orange) and competitor A (red) decrease in stress at 98.6 °F (37 °C) in comparison to 73.4 °F (23 °C).*

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Fig. 41: THERMEO® time depending stress relaxation at 73.4 °F (23 °C).*

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The thousands of times proven **THERMEO®** round is now also available as an **MSI®** (Multi Species Inhibition) version. This makes it the world's first biomimetic milling disc in its class, which naturally reduces the adhesion of biofilms in the oral environment using **pro3dure's** unique **MSI® technology**. It copies a natural process, such as that used by the red alga *Delisea*. A natural substance (lactam) anchored in the dental material interferes with the communication of bacteria, the so-called "quorum sensing". This significantly reduces the formation of biofilms or their adhesion. This effect is not due to antibacterial additives, which may lead to bacterial resistance. In the course of this, improved oral hygiene is achieved and thus less irritation of the soft tissue in accordance with the motto "Inspired by nature, addicted to function".

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Comprehensive documentation and important tips on the use and care of **THERMEO®** splints help you and the patients use **THERMEO®** products for a long and sustainable time.

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“With **THERMEO® | MSI**, processing, function, fit and aesthetics can be perfectly combined!”

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THERMEO® SO milling blanks are used to produce functional aesthetic splints (snap-on splints) with a thermomemory effect and are characterized by their thermoplastic flexibility. This ensures extremely precise adaptation to the tooth situation and exceptional, tension-free wearing comfort for the patient. The self-adjusting **THERMEO® SO** snap-on splints are extremely fracture-resistant and have an optimized translucency in the tooth shades A1, A2 and A3. The industrial polymerization process achieves maximum material homogeneity and thus guarantees outstanding long-term stability. Snap-on splints are temporary restorations for the pre-prosthetic treatment of complex prosthetic cases. **THERMEO® SO** – long-term stable and functional dentures with a high degree of safety, especially for diagnostic occlusal adjustment.

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We would like to take this opportunity to introduce you to a world first from the **pro3dure** development department. With the thermoforming material **ISM-1** for individual sports mouthguards, we have succeeded for the first time in making **MSI® technology** available for EVA-based films. In a study by Prof. Gordon Ramage (Glasgow Caledonian University), it was shown that this new class of material has significantly lower biofilm formation compared to conventional EVA-based mouthguards. This is shown below in **Fig. 42** for *Candida albicans*.

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Colonization and recolonization of *Candida albicans* after long-term incubation (24 h). Conventional EVA mouthguard films (a) and **ISM-1 | MSI** films (b) were examined under a microscope for evaluation.

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Insulating agent for
3D-printed dental models

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The new **printisol** from **pro3dure** is an insulating agent that forms an insulating layer between the 3D-printed model and the cold-curing resin. The cold-curing polymer can be easily lifted off the model and has an excellent fit to the plastic base. Apply with a brush for an optimum insulating effect. For this purpose, high-quality resin models¹ produced using the 3D printing process (e.g. with the Asiga MAX™) should be used.

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“**printisol** is finally a functioning insulation for 3D printed models!”

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- dimensionally stable and accurate fit
- thixotropic, stable and controllable
- low shrinkage
- residue-free combustion

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LCR-1 is a light curing gel which can be used in a variety of applications. These applications include but are not limited to model cast technology, use as block-out gel for bleaching trays and as a repair material of 3D printed cast models. The gel is thixotropic, therefore steady and manageable when modeling. The marginal shrinkage of **LCR-1** leads to exceptional high dimensional stability and accuracy of fit.

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1 l Flasche / 1 l bottle

Art.-Nr. / item-no.: **U1900001**

Water soluble cleaning agent
for the removal of resin residues
on additively manufactured
dental molded parts.

- non explosive
- ideal for ultrasonic systems
- highly efficient
(up to 5 times more than Isopropanole)
- ultrafiltrable
- compatible with **CLD-1**

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The compact 3D printing system for dentistry

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Compact, user-friendly and precise are the attributes of the **Prusa Medical System**, which consists of the **Prusa Medical One** 3D printer and the **CW One** washing and post-curing unit. The Prusa Medical One has a high-resolution UV-LED panel and belongs to the group of image projection systems based on LED-LCD technology. A variety of special solutions, such as build platforms specially adapted to dental requirements and software features to compensate for curing layers, guarantee precise and reproducible results in everyday work. In combination with the all-in-one solution CW One, which allows you to wash and post-cure your components, a complete process chain for the production of dental molded parts can be realized in a minimum of space at an exceptional price level. With this system, you have access to validated processes with **pro3dure's printodent** materials and thus to the world's largest portfolio of dental applications for 3D printing. The indications range from precise dental models with **printodent® GR-13 model**, splints and drill templates with **printodent® GR-10.1 guide | MSI**, denture bases with **printodent® GR 14.2 HI | MSI**, denture teeth, crowns and bridges for temporary long-term care with **printodent® GR-17.1 temporary It**. Prusa Medical One convinces with precise printing results and easy handling.

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- Easy to use
- intuitive menu navigation and optimized workflow
- large build platform (127 x 80 x 150 mm)
- precise components for reproducible accuracy
- unique calibration
- convincing print results
- validated processes

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"With **pro3dure** and **Prusa**, I have all validated parameters directly from the manufacturer. So my workflow is precise and accurate."

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"We are very happy that we have teamed up with **pro3dure**! They provide top-notch biocompatible resins and are excellent business partners."

Josef Prusa
(CEO and Founder Prusa Research)

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Fig. 43: Precision printing platforms are an important part of our Medical One 3D printer. They are manufactured with maximum precision so that they can be freely exchanged without the need to calibrate the printer. This results in a significant reduction of the downtime between prints and even higher efficiency of the entire printing process.

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Fig. 44: The **Prusa Slicer** software package features a large number of special solutions that are specifically tailored to the requirements of the dental sector.

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The new **ASIGA MAX 2™** forms the backbone of your digital process chain with the **printodent®** materials. This 3D printer, specially developed for the digital workflow in the laboratory, is characterized by maximum productivity with minimum space requirements. With a pixel size (x,y) of 62 µm and a 385 nm UV irradiation source, you have access to the complete material portfolio of the **pro3dure printodent® series**. This makes the ASIGA MAX 2™ 3D printer an all-rounder.

Technische Daten technical data		ASIGA MAX 2		Art.-Nr. / item-no.: U2016052
Technologie / technology:	DLP	Wellenlänge / wave length:	385 nm	
Leistung / power:	12 V /120 W	Bauraum / build envelope:	119 x 67 x 75 mm	
Spannungsquelle / power supply:	100-240 V	Gewicht /weight:	19,3 kg	
Lichtquelle / irradiation source:	UV-LED	Laterale Auflösung / lateral resolution:	62 µm	
Maße (B x T x H) / Dimensions (L x W x H):	260 x 380 x 370 mm	Schichtstärke / layer thickness:	25-100	

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
Discover the new **Asiga Ultra™**. The Active Layer Monitoring technology guarantees precise layer thickness control for maximum print accuracy and speed. Experience razor-sharp 4K image resolution over a large build area. With the auto-calibrating UV-LED system and DLP micromirror technology, the Asiga Ultra offers reliability and versatility. Touchless entry allows you to interact effortlessly without contaminating surfaces. The magnetic platform holder and infrared material heater optimize your workflow and enable the use of a wide range of advanced **pro3dure printodent®** materials. The intuitive user interface and Transparent Mode™ make the printing process easier than ever. Rely on chamber heating for optimal process reproducibility.

Technische Daten technical data		ASIGA Ultra (50)		Art.-Nr. / item-no.: U2016053
Technologie / technology:	DLP	Wellenlänge / wave length:	385 nm	
Leistung / power:	500 W	Bauraum / build envelope:	192 x 108 x 147 mm	
Spannungsquelle / power supply:	100-240 V	Gewicht /weight:	32,5 kg	
Lichtquelle / irradiation source:	UV-LED	Laterale Auflösung / lateral resolution:	50 µm	
Maße (B x T x H) / Dimensions (L x W x H):	353 x 441 x 555 mm	Schichtstärke / layer thickness:	25-100	

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Processor: 2.4 GHz multi-core
Memory: 4 GB
Graphic: 256 MB, OpenGL 2.0 or later support
Mouse: 3-button mouse with scroll wheel
Network: Required for printing
Storage: 1 GB free

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Enhance the performance of your **Asiga** system in your dental environment by using the **Asiga Composer software package**. Asiga Composer acts as a link between your CAD-Software and your Asiga printer in order to realize a smooth and integrated digital workflow. The Composer is an intuitive, automatic STL & SLC support and build preparation package. Based on the novel optimal support algorithm, parts are supported automatically only where absolutely required by the part geometry. This makes job setup and post-processing a breeze. With Composer you can include any STL, SLC and Stomp files in the same build. This saves users from having to convert files from one type to another. Projects can be saved as a lightweight "Asiga Composer Build" (.comblld) or complete export into a ZIP file. Composer's automatic part placement tool handles efficient placement of both STL & SLC files based on part con-tours. This maximizes your available build area. You can manage your workflow by preparing more than one print at the same time and switch between platforms which are conveniently tabbed at the top of your screen. Composer is a cross platform desktop software supporting Windows, Macintosh and Linux. It is included free with every system with an unlimited seats license for use throughout your office.

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The **printo[®]dent** family can now also be used for HPS technology **Hybrid PhotoSynthesis™** with the **X1 LUMIA DENTAL**.

Eliminate all tradeoffs between throughput, accuracy, feature resolution and surface finish. The Lumia XI Hi-Speed SLA system is the answer to the challenges of dental 3D printing. It combines unrivaled precision with impressive production capacity without compromising on accuracy, detail resolution and surface quality. The Hybrid PhotoSynthesis (HPS) process offers exceptional precision and accuracy in the production of dental models, splints, surgical guides, and prosthetics at very high throughputs. Contact us about the validated processes!

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The „cleaning device“ **CLD-1** represents another building block for the optimized and controlled manufacturing process of medical devices by generative fabrication. Unlike conventional ultrasonic baths the **CLD-1** generates a vortex of the cleaning agent. As a result, microcracks in the component surfaces are significantly reduced, thus further increasing the part quality in terms of surface finish, appearance and mechanics of your dental objects. A selection of different cleaning programs assists you in perfectly matching the cleaning program to the selected material class and choice of cleaning agent. The cleaning container can be easily removed and changed in order to simplify the manufacturing processes in your lab. Changing the cleaning liquid can be done within seconds in a convenient way.

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The **RLR-1** is a roller/tilting stirring device for homogenization of 3D printing materials with pigment or filler content. Highly filled 3D printing materials like the **printodent® GR-17.1 temporary It** require a thorough mixing before use, handshaking is often insufficient. This can result in significant colour deviations and even build failures may occur. The RLR-1 Mixer is constructed for constant operation in order to keep your **pro3dure** 3D materials ready for use at any time with the same consistency.

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
The high-tech post-curing unit **CD-2** is tuned to the latest generation of 3D printers. It features top-level productivity, the most current LED irradiation technology, sophisticated sensors and a wide range of smart features. Through the touch panel you can fully control and monitor all relevant irradiation parameters. Sensors measure the protective gas flow and inform the user about the status of the light source. This avoids inhibition layers on generated parts. As a result, the biocompatibility and the optimal mechanical properties of your 3D printed parts are guaranteed. In a profound and unprecedented form the CD-2 combines compactness and efficiency on the highest level of technology, along with an excellent price/performance ratio.

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printo[®]
dent

		GR-10 guide	GR-10 guide MSI	GR-10.1 guide MSI	GR-11 tray	GR-11.1 tray	GR-12 cast	GR-13 model	GR-13.1 model	GR-13.2 model	GR-13.3 model release	GR-13.4 model	GR-14.2 denture HI	GR-14.3 denture HI MSI	GR-15 gingiva	GR-16 Xray	GR-17 temporary	GR-17.1 temporary It	GR-18.2 IB	GR-19.1 OA	GR-19.1 OA MSI	GR-20 MJF	GR-21.1 Try-In	GR-22 flex	GR-23 mouthguard
ASIGA MAX		✓	○	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	○	✓	✓	✓	✓	✓
ASIGA MAX UV / MAX 2		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
ASIGA PRO 2		✓	○	○	✓	○	✓	✓	✓	○	✓	○	○	✓	✓	✓	○	○	○	✓	○	✓	✓	○	○
ASIGA Ultra		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
ASIGA PRO 4K		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
SprintRay Pro 95		○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○
MiiCraft Ultra		✓	✓	✓	✓	✓	✓	✓	✓	○	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	○	○
MiiCraft Prime		✓	✓	✓	✓	✓	✓	✓	✓	○	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	○	○
MiiCraft Profession 4K		✓	✓	✓	✓	✓	✓	✓	✓	○	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	○	○
MiiCraft Alpha		✓	✓	✓	○	✓	✓	✓	✓	○	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	○	○
Organical X1S		✓	○	○	✓	○	✓	✓	✓	○	○	✓	○	○	✓	✓	○	○	○	✓	○	✓	○	○	○
Organical X1N		✓	✓	✓	✓	✓	✓	✓	✓	○	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	○	○
Organical 3D Print X1P-9K		✓	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○
phrozen Sonic XL 4K		✓	○	○	✓	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○

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printc
dent

		GR-10 guide	GR-10 guide MSI	GR-10.1 guide MSI	GR-11 tray	GR-11.1 tray	GR-12 cast	GR-13 model	GR-13.1 model	GR-13.2 model	GR-13.3 model release	GR-13.4 model	GR-14.2 denture HI	GR-14.3 denture HI MSI	GR-15 gingiva	GR-16 Xray	GR-17 temporary	GR-17.1 temporary	GR-18.2 IB	GR-19.1 OA	GR-19.1 OA MSI	GR-20 MJF	GR-21.1 Try-In	GR-22 flex	GR-23 mouthguard
Flashforge		✓	✓	✓	○	✓	✓	✓	✓	✓	○	✓	✓	✓	✓	✓	✓	✓	✓	○	✓	✓	○		
DENTAL MAKERS		✓	○	○	○	○	○	○	✓	○	○	✓	○	○	○	✓	○	○	○	○	✓	○	○		
RAYSHAPE		✓	○	○	○	○	○	○	○	○	○	✓	○	○	○	○	✓	○	✓	○	○	○	○		
UniFormation		✓	○	○	✓	○	○	○	○	○	○	✓	✓	○	○	✓	○	○	✓	○	✓	○	○		
ZYLO 3D		✓	○	○	○	○	✓	○	○	○	○	○	○	○	○	○	✓	○	✓	✓	○	○	✓	○	

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printc
dent

		printo dent																							
		GR-10 guide	GR-10 guide MSI	GR-10.1 guide MSI	GR-11 tray	GR-11.1 tray	GR-12 cast	GR-13 model	GR-13.1 model	GR-13.2 model	GR-13.3 model release	GR-13.4 model	GR-14.2 denture HI	GR-14.3 denture HI MSI	GR-15 denture HI MSI	GR-16 Xray	GR-17 temporary	GR-17.1 temporary	GR-18.2 IB	GR-19.1 OA	GR-19.1 OA MSI	GR-20 MJF	GR-21.1 Try-In	GR-22 flex	GR-23 mouthguard
pro3dure		✓	✓	✓	✓	✓	✓	✓	✓	✓	○	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	○	
SprintRay		○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	
Kulzer		✓	✓	✓	✓	✓	✓	✓	✓	○	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	○		
Prusa		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		
rapidshape		✓	✓	✓	○	✓	✓	✓	✓	✓	○	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	○		
Meccatronicore		✓	✓	✓	○	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	○		

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printc[®]
dent

The largest portfolio of 3D printing
resins for dental applications by ...

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printo
dent

printo dent		GR-10 guide GR-10 guide MSI GR-10.1 guide MSI GR-11 tray GR-11.1 tray GR-12 cast GR-13 model GR-13.1 model GR-13.2 model GR-13.3 model release GR-13.4 model GR-14.2 denture HI GR-14.2 denture HI MSI GR-14.3 denture HI MSI GR-15 gingiva GR-16 Xray GR-17 temporary GR-17.1 temporary It GR-18.2 IB GR-19.1 OA GR-19.1 OA MSI GR-20 MJF GR-21.1 Try-In GR-22 flex GR-23 mouthguard																							
pro3dure CD-1		✓	✓	✓	✓	✓	✓	✓	✓	✓	○	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	○	
pro3dure CD-2		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
NK-OPTIK Otoflash G171		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
ASIGA Flash		○	○	✓	○	✓	✓	✓	✓	✓	○	○	○	○	○	○	○	○	○	○	○	○	○	○	
ASIGA Cure 1.5 I		✓	○	○	○	✓	✓	✓	✓	✓	✓	✓	○	✓	✓	✓	✓	○	✓	✓	✓	✓	✓	○	
RAYSHAPE ShapeCure D		✓	○	○	○	○	○	○	○	○	○	✓	○	○	○	○	○	✓	○	✓	✓	○	○	○	
SprintRay ProCure		○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	
SprintRay ProCure 2		○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	
SprintRay NanoCure		○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	
SISMA EVERES CURO		✓	✓	✓	✓	✓	✓	✓	✓	✓	○	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	○	
Prusa MEDICAL CW ONE		✓	✓	✓	✓	✓	✓	✓	✓	✓	○	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	○	
WICKED ENGINEERING CUREbox		✓	✓	✓	✓	✓	✓	✓	✓	✓	○	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	○	

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GR-22 flex

GR-23 mouthguard

rapidshape RS Cure 	✓	✓	✓	○	✓	✓	✓	✓	✓	✓	○	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	○	
Nobil Metal Dentafarm Photopol 	✓	✓	✓	○	✓	✓	✓	✓	✓	○	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	○	
Meccatroncore BB Cure N 	✓	✓	✓	○	✓	✓	✓	✓	✓	✓	○	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
Dreve PCU LED N2 	✓	○	○	○	✓	✓	✓	○	○	○	○	✓	○	○	○	○	○	○	✓	○	✓	✓	○	○	○
Kulzer cara Print LEDCure 	✓	✓	✓	✓	✓	✓	✓	✓	✓	○	✓	✓	✓	○	✓	✓	✓	○	✓	✓	○	✓	✓	○	
SHINING 3D FabCure 	✓	✓	✓	✓	✓	✓	✓	✓	✓	○	✓	✓	✓	✓	✓	✓	○	✓	✓	✓	○	✓	✓	○	
SHINING 3D FabCure N2 	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	✓	○	✓	✓	○	○	○	

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All information contained in this brochure has been compiled with the utmost care and reflects the status as of the date of publication. We reserve the right to make changes as part of ongoing technical development. The information provided in the current instructions for use and official product specifications shall always take precedence.

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pro3dure medical LLC