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

**SysTherm® PU** is a 2-component polymer for the production of earmolds with thermoeffect. The cured **SysTherm® PU** material is characterized by an initial Shore A hardness of 85 at room temperature (23 °C) which changes to 65 Shore A at body temperature (37 °C). Accordingly, **SysTherm® PU** earmolds can be worked on at room temperature. Moreover, the end customer can precisely position a **SysTherm®** earmold in their ear canal. At body temperature the earpiece becomes soft and flexible, generating outstanding wearing comfort in combination with a perfect sealing during jaw movements. **SysTherm®**'s excellent recovery ability ensures that the earmold regains its original shape after deformation. When using a digital workflow the **SysTherm® PU** material can be processed with 3D printed castforms (pro3dure's GR-12.2) and the **SysTherm® Separator** or solely with pro3dure's **SysTherm® Cast**. Additionally, the usage of the material in the traditional PnP process with a silicone negative form (pro3dure's FS-2) is also possible. Due to its chemical base **SysTherm®** is moisture sensitive and should be stored dry at room temperature (23 °C).

- 1) Open the **SysTherm®** secondary packaging (fig. 1).
- 2) Twist off the cartridge top and attach the mixing canula (fig. 2).
- 3) Insert the cartridge into the **SysTherm®**-dispenser (fig. 3).
- 4) Before dispensing into a mold please pump a pearl-sized amount of **SysTherm® PU** through the mixing tip in order to ensure a perfect mixture of the material (fig. 4).
- 5) Depending on the chosen production process please proceed to one of the corresponding sections below.
- 6) For storing of an opened **SysTherm® PU** cartridge please leave the mixing canula on the cartridge (fig. 6).

### A)


Production process based on 3D-printed molds (GR 12.2) with cast separator (**SysTherm® Separator**):

- 1) Printing cast shells out of GR 12.2 (Make sure air-vents are the same height as filling inlet, otherwise PU material will leak out during curing) (fig. 7).
- 2) After printing thoroughly clean the parts using IPA and then fresh IPA, for 5 min. each (fig. 8).
- 3) Dry the parts with compressed air (fig. 9).
- 4) Dry the parts in oven at 40 °C for at least 30 min. to eliminate residual IPA. Parts must be odor free.
- 5) Cure the parts with nitrogen for 2 min. in **CD-2**.
- 6) Dip the casts into separator liquid – **SysTherm® Separator** (fig. 10).
- 7) Thoroughly dry the cast with compressed air (fig. 11).
- 8) Dry in oven at 60 °C for 120 min. (fig. 14).
- 9) Thoroughly dry the cast with compressed air.  
Note: the separator will stay as a thin liquid film on the surface (fig. 11).
- 10) Fill in **SysTherm® PU** carefully and avoid trapped air. Fill the inlet as a reservoir for possible polymerization shrinkage (fig. 13).
- 11) Put the filled casts into a dry pressure pot – make sure **SysTherm®** will not flow out.
- 12) Cure **SysTherm®** at 60 °C for 2 hours. Keep heating up time into account (fig. 14).
- 13) Cool down the filled parts and carefully crack the shells. Use eye protection and appropriate tools (fig. 15).
- 14) Peel off cracked cast parts.
- 15) Clean and finish **SysTherm®** parts.

Indications for use:  
silicone earmold material

### Technical data::

- Mixing volume:  
50 ml cartridge
- Mixing ratio:  
4:1
- Setting time:  
120 min. at  
60 °C / 140 °F  
in pressure pot
- Final hardness:  
Shore A: 23°:85 ±5  
Shore A: 37°:70 ±5
- Application:  
at 23 ± 2 °C,  
50 ± 5 % rel.  
humidity
- Storage:

### Ordering information:

**SysTherm® Polymer**  
2 x 50 g  
10 mixing canulas  
Art.-Nr.: A4020124

**SysTherm® Cast**  
1000 g  
Art.-Nr.: A4020125


**SysTherm® SPT-1**  
500 g  
Art.-Nr.: A4020126

**SysTherm® mixing canulas**  
25 pieces  
Art.-Nr.: A4020127

**SysTherm® 4:1 dispenser**  
1 piece  
Art.-Nr.: A4020128

CE  
0481

These data result from measurements of a representative sample, which were determined within the scope of our quality assurance.

### B)


Production Process based on 3D-printed molds (**SysTherm®**):

- 1) Printing cast shells out of **SysTherm® Cast**. Make sure the air-vents are the same height as filling inlet, otherwise PU material will leak out during curing (fig. 7).
- 2) After printing thoroughly dry the parts with compressed air to remove residual liquid resin (fig. 9).
- 3) Cure the parts with nitrogen for 2 min. in **CD-2**.
- 4) Do NOT use any kind of cast separators.
- 5) Dry the casts in oven at 40 °C for 60 min. to remove any moisture (fig. 12).
- 6) Fill in **SysTherm® PU** carefully and avoid trapped air. Fill the inlet as a reservoir for possible polymerization shrinkage (fig. 13).
- 7) Put the filled casts into a dry pressure pot. Make sure **SysTherm®** will not flow out.
- 8) Cure **SysTherm®** at 60 °C for 2 hours. Keep heating up time into account (fig. 14).
- 9) Place the casts into a beaker with hot water (60-70 °C) and wait 30 min. Cast material will swell (fig. 16).
- 10) Peel off swelled cast material and remove parts (fig. 15).
- 11) Clean and finish **SysTherm®** parts.

### C)

Production process based on PNP method using form silicones :

- 1) Preparation of negative mould out of FS-2 silicone (see instruction for use).
- 2) If necessary cleaning of mould and drying (60 °C, 120 min.) to avoid moisture (fig. 12).
- 3) Fill in **SysTherm® PU** carefully and avoid trapped air (fig. 18).
- 4) Put the filled casts into a dry pressure pot – make sure **SysTherm®** will not flow out.
- 5) Cure **SysTherm®** at 60 °C for 2 hours. Keep heating up time into account (fig. 14).
- 6) Clean and finish **SysTherm®** parts.

Temperature limit

CE mark

Consult instructions for use

Use by date

Batch code

Catalogue number

Date of manufacture

idea to product.

Made in Germany  
pro3dure medical GmbH  
Am Burgberg 13  
D - 58642 Iserlohn  
Tel. +49 (0)2374 920050-0  
Fax +49 (0)2374 920050-50  
info@pro3dure.com  
www.pro3dure.com