

Earmould Silicone ES-1.1/25

Instructions

pro3dure
medical


pro3dure medical GmbH

Am Burgberg 13
D - 58642 Iserlohn
Tel. +49 (0)2374 920050-0
Fax +49 (0)2374 920050-50
info@pro3dure.com
www.pro3dure.com

Indications for use:
silicone earmould
material

Technical data:







- Mixing volume:
48 ml cartridge
- Mixing ratio: 1:1
- Color code:
Base: depends on
chosen material
Catalyst: white
- Setting time:
60 min. at
50 °C / 122 °F
in pressure pot
- Final hardness:
25 ± 5 Shore A
- Application:
At 73 ± 4 °F,
50 ± 5 % rel.
humidity
- Storage:



**Ordering
information:**

cartridges of 8 x 48 ml
clear-transparent
item no.: A0030800
reddish-transparent
item no.: A0030203
black-opaque
item no.: A0030376
white-opaque
item no.: A0030100
neon-green-opaque
item no.: A0030530
neon-yellow-opaque
item no.: A0030729
neon-pink-opaque
item no.: A0030903
neon-orange-opaque
item no.: A0031327
dark-blue-opaque
item no.: A0030428
light-blue-opaque
item no.: A0031128
red-opaque
item no.: A0030627
violet
item no.: A0030680
smoke grey
item no.: A0030690
Injector DS
50 ml / Stück
item no.: A40000
5,4 mm
Clear Mixing Tip
100 Stück
item no.: A40002

These data result from
measurements of a
representative sample,
which were determined
within the scope of our
quality assurance.



The **ES-1.1/25** is a permanently elastic addition-vulcanizing silicone material for the generation of soft silicone and hearing protection earmoulds.

1. Mixing and dispensing

Place the cartridge (**fig 1**) into the dispensing gun. Unscrew cap. In order to ensure a proper flow and mixing ratio from both orifices (**fig 2**) extrude carefully a small amount of the material. Insert the static mixer into the guided grooves of cartridge. Rotate the mixing cannula anti-clockwise to a stop (**fig 3**). If necessary, mount an Intra-tip. For reducing forces during dispensing, static mixers of a larger diameter (9mm) can be used alternatively. Dispenser is now ready to mix the 2 component earmould silicone in any amount needed.

2. Application (PnP method)

In combination with the **ES-1.1/25** material all commercially available materials for negative forms like plasters, gels or 3D-printing resins can be used for the manufacturing of the negative form (**fig 4**). If necessary, coat the negative form with a separation liquid (e.g. alginate based for plaster). 3D-printed cast forms have to be cleaned carefully before use in order to avoid inhibition of the silicone reaction. Inject the **ES-1.1/25** material slowly and bubblefree into the generated negative form (**fig.5**). In case of narrow structures (e.g. auditory canals), we recommend using a thin mixing canula. Put the negative mold with the injected material into a pressure pot for about 60 min at 122 °F (50 °C).

3. Final surface treatment and lacquering

The silicone earmould can now be shaped by using special cutters and grinding sleeves. Roughen the surface with corundum paper (grit 180) before varnishing. Clean the surface and coat it with a silicone lacquer according to the special instructions of the lacquer (**fig 6**).

4. Important working hints

- In case of the generative fabricated castforms a careful preparation as well as a thorough cleaning will ensure better and more reliable results.
- Ensure that all traces of residual resin is removed from the negative form.
- The cured earmould material **ES-1.1/25** is chemically inert.
Avoid spots on clothing.
- We recommend standard gloves made of nitrile or polyethylene.

Safety advice:

pro3dure medical is not liable for any damage caused by improper application of the earmould material.

Despite of the fact that the high biocompatibility of silicone materials is proven over decades in exceptional cases undesired reactions of the immune system like allergies, irritations can't be absolutely excluded. In case of doubt, we recommend to contact your consulting physician and make an allergy test before the application of the material.

For use by trained specialists.

ES-1.1

idea to product.

CE
0481