

## SECTION 1: Identification of the substance/mixture and of the company/undertaking

### 1.1. Product identifier

**RPR-1**

UFI

QFJF-HUHR-R119-7VN8

### 1.2. Relevant identified uses of the substance or mixture and uses advised against

**Relevant identified uses:** blocking-out, modelling

**Uses advised against:** Private households (= general public). Observe instructions for use.

### 1.3. Details of the supplier of the safety data sheet

#### Manufacturer

pro3dure medical GmbH

Am Burgberg 13  
D 58642 Iserlohn

**Telephone** +49 (0)2374 920050-10

**Telefax:**

#### Supplier

pro3dure medical GmbH

Am Burgberg 13  
D 58642 Iserlohn

**Telephone** +49 (0)2374 920050-10

**Telefax:**

#### Information contact

pro3dure medical GmbH

**Information telephone** +49 (0)2374 920050-10

**Information telefax**

**E-mail (competent person)** info@pro3dure.com

**Website** www.pro3dure.com

### 1.4. Emergency telephone number

pro3dure medical GmbH  
This number is serviced during office hours.

**Telephone** +49 (0)2374 920050-10

## SECTION 2: Hazards identification

### 2.1. Classification of the substance or mixture

Regulation (EC) No 1272/2008:  
Skin Sens. 1, H317; Repr. 1B, H360

### 2.2. Label elements

Classification according to Regulation (EC) No 1272/2008 [CLP]

Hazard pictograms



GHS07, GHS08  
Danger

**Signal word:**

**Hazard statements:**

H317 May cause an allergic skin reaction.  
H360 May damage fertility or the unborn child.

**Precautionary statements:**

P201 Obtain special instructions before use.  
P202 Do not handle until all safety precautions have been read and understood.  
P261 Avoid breathing dust/fume/gas/mist/vapours/spray.  
P263 Avoid contact during pregnancy and while nursing.  
P264.1 Wash hands thoroughly after handling.  
P272 Contaminated work clothing should not be allowed out of the workplace.  
P280 Wear protective gloves/protective clothing/eye protection/face protection.  
P302+352.1 IF ON SKIN: Wash with plenty of soap and water.  
P305+351+338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.  
P308+313 IF exposed or concerned: Get medical advice/attention.  
P333+313 If skin irritation or rash occurs: Get medical advice/attention.  
P362+364 Take off contaminated clothing and wash it before reuse.  
P403+233 Store in a well-ventilated place. Keep container tightly closed.  
P405 Store locked up.  
P501 Dispose of contents/container in accordance with local and national regulations.

**Special labelling of particular preparations:**

Hazardous component(s) for labelling: Ethyl 4-dimethylaminobenzoate, Diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide, 7,7,9(or 7,9,9)-trimethyl-4,13-dioxo-3,14-dioxo-5,12-diazahexadecane-1,16-diyl bismethacrylate

## 2.3. Other hazards

The substances in the mixture do not meet the PBT/vPvB criteria according to REACH, annex XIII.  
No endocrine-disrupting properties known (see section 12)  
Annex XVII, entry: 75, 78 not relevant  
Contains: SVHC substance.  
CA Proposition 65: Substance(s) not listed.

## SECTION 3: Composition / information on ingredients

### 3.1. Substances

not applicable

### 3.2. Mixtures

Mixture with, among others, the following ingredients and other non-hazardous admixtures

**Composition/information on ingredients**

| Substance:                                                                                     | CAS-No.:   | REACH-no.:            | Concentration: | Classification: EC 1272/2008 (CLP):       | M, ATE, Note                                                                                                  |
|------------------------------------------------------------------------------------------------|------------|-----------------------|----------------|-------------------------------------------|---------------------------------------------------------------------------------------------------------------|
| Diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide                                                | 75980-60-8 | 01-2119972295-29-xxxx | <02            | Skin Sens. 1, H317; Repr. 1B, H360Fd      | M = 0 ATE (dermal) = > 2000 mg/kg bw ATE (oral) = > 5000 mg/kg bw ATE (inhalativ) = -                         |
| 7,7,9(or 7,9,9)-trimethyl-4,13-dioxo-3,14-dioxo-5,12-diazahexadecane-1,16-diyl bismethacrylate | 72869-86-4 | 01-2120751202-68-XXXX | < 80           | Skin Sens. 1, H317                        | M = 0 ATE (dermal) = >2000 mg/kg bw ATE (oral) = >5000 mg/kg bw ATE (inhalativ) = - Classification acc. ATP21 |
| Ethyl 4-dimethylaminobenzoate                                                                  | 10287-53-3 |                       | < 0,5          | Repr. 1B, H360FD; Aquatic Chronic 2, H411 | M = 0 ATE (dermal) = - ATE (oral) = 2000 mg/kg ATE (inhalativ) = -                                            |

(Full text of H- and EUH-statements: see section 16.)

## SECTION 4: First aid measures

#### 4.1. Description of first aid measures

|                                |                                                                                                                                                                                        |
|--------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>General information:</b>    | In case of accident or unwellness, seek medical advice immediately Remove contaminated, saturated clothing immediately. When in doubt or if symptoms are observed, get medical advice. |
| <b>In case of inhalation:</b>  | Provide fresh air. Seek medical attention if problems persist.                                                                                                                         |
| <b>Following skin contact:</b> | After contact with skin, wash immediately with plenty of water and soap. In case of skin irritation, consult a physician.                                                              |
| <b>After eye contact:</b>      | IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Seek medical attention if problems persist.           |
| <b>After ingestion:</b>        | Rinse mouth immediately and drink plenty of water. No administration in cases of unconsciousness or cramps. Do NOT induce vomiting. Call a physician immediately.                      |

#### 4.2. Most important symptoms and effects, both acute and delayed

Adverse human health effects and symptoms : May cause an allergic skin reaction. erythema (redness). Causes skin irritation.

#### 4.3. Indication of any immediate medical attention and special treatment needed

First Aid, decontamination, treatment of symptoms. May cause sensitisation especially in sensitive humans. May cause sensitisation especially in sensitive humans. Provide fresh air.

---

### SECTION 5: Firefighting measures

#### 5.1. Extinguishing media

|                                       |                                                                          |
|---------------------------------------|--------------------------------------------------------------------------|
| <b>Suitable extinguishing media</b>   | ABC-powder. BC-powder. alcohol resistant foam. Carbon dioxide. Nitrogen. |
| <b>Unsuitable extinguishing media</b> | Water spray jet. Full water jet. Excess water.                           |

#### 5.2. Special hazards arising from the substance or mixture

In case of fire may be liberated: Carbon dioxide (CO<sub>2</sub>). Carbon monoxide.

#### 5.3. Advice for firefighters

##### General information

Move undamaged containers from immediate hazard area if it can be done safely. Use water spray jet to protect personnel and to cool endangered containers. Do not allow water used to extinguish fire to enter drains or waterways. Co-ordinate fire-fighting measures to the fire surroundings.

##### Special protective equipment for fire-fighters:

In case of fire: Wear self-contained breathing apparatus. Wear chemical resistant suit.

---

### SECTION 6: Accidental release measures

#### 6.1. Personal precautions, protective equipment and emergency procedures

Provide adequate ventilation. Remove all sources of ignition. Remove persons to safety. Wear personal protection equipment. Special danger of slipping by leaking/spilling product. See protective measures under point 7 and 8.

#### 6.2. Environmental precautions

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil.

#### 6.3. Methods and material for containment and cleaning up

Absorb with liquid-binding material (e.g. sand, diatomaceous earth, acid- or universal binding agents). Treat the recovered material as prescribed in the section on waste disposal.

## 6.4. Reference to other sections

Safe handling: see section 7 Personal protection equipment: see section 8 Disposal: see section 13

## SECTION 7: Handling and storage

### 7.1. Precautions for safe handling

#### Advices on safe handling

Keep container tightly closed. When using do not eat, drink, smoke, sniff. Use only in well-ventilated areas. All work processes must always be designed so that the following is as low as possible: Inhalation Avoid contact with skin, eyes and clothes.

#### Precautions against fire and explosion:

Always close containers tightly after the removal of product. Use only in well-ventilated areas. Usual measures for fire prevention.

### 7.2. Conditions for safe storage, including any incompatibilities

#### Requirements for storage rooms and vessels

Keep container tightly closed in a cool, well-ventilated place. Protect against: UV-radiation/sunlight Only use the material in places where open light, fire and other flammable sources can be kept away.

#### Hints on joint storage

Do not store together with: Oxidizing agent. Organic peroxides. Keep away from food, drink and animal feedingstuffs.

### 7.3. Specific end use(s)

Observe instructions for use.

## SECTION 8: Exposure controls/personal protection

### 8.1. Control parameters

#### occupational exposure limit value

| Substance: | CAS-No.: |  | Source: | Occupational exposure limit value:[ppm] | Occupational exposure limit value:[mg/m <sup>3</sup> ] | Limitation of exposure peaks: | Remark: |
|------------|----------|--|---------|-----------------------------------------|--------------------------------------------------------|-------------------------------|---------|
|------------|----------|--|---------|-----------------------------------------|--------------------------------------------------------|-------------------------------|---------|

#### Substance with a common (EC) occupational exposure limit value.

| Substance: | CAS-No.: |  | Source: | Occupational exposure limit value:[ppm] | Occupational exposure limit value:[mg/m <sup>3</sup> ] | Limitation of exposure peaks: | Remark: |
|------------|----------|--|---------|-----------------------------------------|--------------------------------------------------------|-------------------------------|---------|
|------------|----------|--|---------|-----------------------------------------|--------------------------------------------------------|-------------------------------|---------|

#### DNEL-/PNEC-values

DNEL value

| Substance:                                                                                      | CAS-No.:   | DNEL/DMEL                                                                                                                                                                                                                                                                                                          |
|-------------------------------------------------------------------------------------------------|------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 7,7,9(or 7,9,9)-trimethyl-4,13-dioxo-3,14-dioxo-5,12-diazaheptadecane-1,16-diyl bismethacrylate | 72869-86-4 | worker inhalative long-term, systemic 3,3 mg/m <sup>3</sup><br>worker dermal long-term, systemic 1,3 mg/kg bw/day<br>population inhalative long-term, systemic 0,6 mg/m <sup>3</sup><br>population dermal long-term, systemic 0,7 mg/kg bw/day<br>population oral long-term, systemic 0,3 mg/kg bw/day             |
| Ethyl 4-dimethylaminobenzoate                                                                   | 10287-53-3 | worker inhalative long-term, systemic 1,65 mg/m <sup>3</sup><br>worker dermal long-term, systemic 2,3 mg/kg bw/day<br>population inhalative long-term, systemic 0,29 mg/m <sup>3</sup><br>population dermal long-term, systemic 0,83 mg/kg bw/day<br>population oral long-term, systemic 0,08 mg/kg bw/day         |
| 3,6,9-trioxaundecamethylene dimethacrylate                                                      |            | worker inhalative long-term, systemic 48,5 mg/m <sup>3</sup><br>worker dermal long-term, systemic 13,9 mg/kg bw/day<br>population inhalative long-term, systemic 14,5 mg/m <sup>3</sup><br>population dermal long-term, systemic 8,33 mg/kg bw/day<br>population oral long-term, systemic 8,33 mg/kg bw/day        |
| Diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide                                                 | 75980-60-8 | worker inhalative long-term, systemic 0,822 mg/m <sup>3</sup><br>worker dermal long-term, systemic 0,233 mg/kg bw/day<br>population inhalative long-term, systemic 0,145 mg/m <sup>3</sup><br>population dermal long-term, systemic 0,0833 mg/kg bw/day<br>population oral long-term, systemic 0,0833 mg/kg bw/day |

### PNEC Value

| Substance:                                                                                      | CAS-No.:   | PNEC                                                                                                                                                                                                                                                                                    |
|-------------------------------------------------------------------------------------------------|------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 7,7,9(or 7,9,9)-trimethyl-4,13-dioxo-3,14-dioxo-5,12-diazaheptadecane-1,16-diyl bismethacrylate | 72869-86-4 | aquatic, freshwater 0,01 µg/l<br>aquatic, marine water 0,001 µg/l<br>sewage treatment plant 3,61 µg/l<br>sediment, freshwater 4,56 µg/kg dw<br>sediment, marine water 0,46 µg/kg dw<br>soil 0,91 mg/kg dw                                                                               |
| Ethyl 4-dimethylaminobenzoate                                                                   | 10287-53-3 | aquatic, freshwater 2 µg/l<br>intermittent release, freshwater 1,9 µg/l<br>aquatic, marine water 0 µg/l<br>intermittent release, marine water 2 µg/l<br>sewage treatment plant 1460 µg/l<br>sediment, freshwater 113 µg/kg dw<br>sediment, marine water 11 µg/kg dw<br>soil 21 µg/kg dw |
| Diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide                                                 | 75980-60-8 | aquatic, freshwater 1,4 µg/l<br>aquatic, marine water 0,14 µg/l<br>sediment, freshwater 115 µg/kg dw<br>sediment, marine water 11,5 µg/kg dw<br>soil 22,2 µg/kg dw                                                                                                                      |

### Additional information

Occupational exposure limit values: No data available, To follow: Toxicological information

## 8.2. Exposure controls

### Occupational exposure controls:

Provide adequate ventilation as well as local exhaust at critical locations. Technical measures and the application of suitable work processes have priority over personal protection equipment.

### General protection and hygiene measures:

When using do not eat, drink, smoke, sniff. Wash hands before breaks and after work. Apply skin care products after work. Wash contaminated clothing prior to re-use. Wear suitable protective clothing, gloves and eye/face protection. Do not breathe dust.

### Personal protection equipment

Only wear fitting, comfortable and clean protective clothing. Further remarks: Minimum standard for preventive measures while handling with working materials are specified in the TRGS 500.

### Respiratory protection

In case of inadequate ventilation wear respiratory protection. Wear breathing apparatus if exposed to vapours/dusts/aerosols. Suitable respiratory protection apparatus: Combination filter device (DIN EN 141). Dust formation : particulates filter device (DIN EN 143). With correct and proper use, and under normal conditions, breathing protection is not required.

### Hand protection

Tested protective gloves are to be worn: DIN-/EN-Norms: EN ISO 374 Suitable material: Butyl caoutchouc (butyl rubber) , NBR (Nitrile rubber). Thickness of the glove material, Breakthrough times and swelling properties of the material must be taken into consideration.

#### **Eyeface protection**

Wear eye/face protection.

#### **Body protection:**

Personal protection equipment. lab coat. Only wear fitting, comfortable and clean protective clothing.  
Barrier creams are not substitutes for body protection.

#### **Environmental exposure controls**

refer to chapter 7. No further action is necessary.

#### **Consumer exposure controls**

refer to chapter 7. No further action is necessary.

#### **Exposure Scenario:**

Skin contact Inhalation

## **SECTION 9: Physical and chemical properties**

### **9.1. Information on basic physical and chemical properties**

#### **Appearance**

#### **Physical state:**

Liquid

#### **Colour:**

According to the product specification

#### **Odour:**

-

#### **Odour threshold:**

-

#### **Safety relevant basis data**

|                                                 | parameter | Value      | unit | Remark                                 |
|-------------------------------------------------|-----------|------------|------|----------------------------------------|
| <b>Melting point/freezing point:</b>            |           |            |      | not determined                         |
| <b>Initial boiling point and boiling range:</b> |           |            |      | not determined                         |
| <b>Flammability:</b>                            |           |            |      | not determined                         |
| <b>lower flammability or explosive limits:</b>  |           |            |      | not determined                         |
| <b>Upper flammability or explosive limits:</b>  |           |            |      | not determined                         |
| <b>Flash point:</b>                             | >         | 140        |      | Lowest value (based on the components) |
| <b>Ignition temperature:</b>                    |           |            |      | not determined                         |
| <b>Decomposition temperature:</b>               |           |            |      | not determined                         |
| <b>pH:</b>                                      |           |            |      | not determined                         |
| <b>Kinematic viscosity:</b>                     |           | 7000-12000 | Pa*s | 23 °C, 10/s                            |
| <b>Water solubility (g/L):</b>                  |           |            |      |                                        |
| <b>Partition coefficient: n-octanol/water:</b>  | Log KOW   | 1.39-3.39  |      | Based on data of the components        |
| <b>Vapour pressure:</b>                         |           |            |      | not determined                         |
| <b>Density:</b>                                 |           | 1,09       | g/mL | not determined                         |
| <b>Relative density:</b>                        |           |            |      | not determined                         |
| <b>Particle properties:</b>                     |           |            |      | not determined                         |

### **9.2. Other information**

Self-accelerating polymerization temperature (SAPT): not determined

## SECTION 10: Stability and reactivity

### 10.1. Reactivity

Can polymerise exothermically if heated, exposed to air, sunlight or by addition of free radical initiators. No hazardous reaction when handled and stored according to provisions.

### 10.2. Chemical stability

The substance is chemically stable under recommended conditions of storage, use and temperature. Self-accelerating polymerization temperature (SAPT): Carriage in accordance with 2.2.41.1.21, not applicable

### 10.3. Possibility of hazardous reactions

Polymerisation.

### 10.4. Conditions to avoid

Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. UV-radiation/sunlight. Do not expose to temperatures exceeding 60 °C/140 °F.

### 10.5. Incompatible materials

Polymerization with heat evolution may occur in the presence of radical forming substances (e.g. peroxides), reducing substances, and/or heavy metal ions.

### 10.6. Hazardous decomposition products

In case of fire may be liberated: Carbon dioxide (CO<sub>2</sub>). Carbon monoxide.

## SECTION 11: Toxicological information

### 11.1. Information on toxicological effects

There are no data available on the preparation/mixture itself. Information comes from reference works, the literature or safety data sheets of the components. Data apply to the component with the highest toxicological risk.

**M-factor:** - **Acute toxicity (dermal):** > 2000 mg/kg  
**Acute toxicity (oral):** > 2000 mg/kg **Acute toxicity (inhalative):** -

#### Acute toxicity

| Substance:                                                                                     | CAS-No.:   | Toxicological information                                                                                                                                                             |
|------------------------------------------------------------------------------------------------|------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide                                                | 75980-60-8 | LC50 inhalation (Rat) 2000 mg/ kg bw<br>LD50 oral (rat) > 5000 mg/kg bw<br>LD50 dermal (rat) > 2000 mg/kg bw<br>NOAEL (rat) 50 mg/kg bw/Tag<br>NOAEL Repr.-Tox. (Rat) 60 mg/kg bw/Tag |
| 3,6,9-trioxaundecamethylene dimethacrylate                                                     |            | LD50 dermal (rat) > 3000 mg/kg<br>LD50 oral (rat) > 5000 mg/kg                                                                                                                        |
| Ethyl 4-dimethylaminobenzoate                                                                  | 10287-53-3 | LD50 oral (rat) > 2000 mg/kg<br>NOAEL (rat) 74 mg/kg bw/Tag<br>NOAEL Repr.-Tox. (Rat) 50 mg/kg bw/Tag                                                                                 |
| 7,7,9(or 7,9,9)-trimethyl-4,13-dioxo-3,14-dioxo-5,12-diazahexadecane-1,16-diyl bismethacrylate | 72869-86-4 | LD50 oral (rat) > 5000 mg/kg<br>LD50 dermal (rat) > 2000 mg/kg<br>NOAEL Repr.-Tox. (Rat) 1000 mg/kg/d<br>NOAEL STOT-RE (rat) 100 mg/kg/d                                              |

#### Skin corrosion/irritation:

Frequently or prolonged contact with skin may cause dermal irritation. May cause an allergic skin reaction.

#### Serious eye damage/irritation:

Irritant.

#### Respiratory or skin sensitisation:

May cause sensitization by skin contact. May cause sensitisation especially in sensitive humans. (based on the components)

### CMR effects (carcinogenicity, mutagenicity and toxicity for reproduction)

#### Carcinogenicity:

There are no data available on the preparation/mixture itself. Based on available data, the classification criteria are not met. (based on the components)

#### Germ cell mutagenicity:

There are no data available on the preparation/mixture itself. Based on available data, the classification criteria are not met. (based on the components)

#### Reproductive toxicity:

There are no data available on the mixture itself. Based on data of the components: NOAEL: 50-1000 mg/kg bw/day oral Rat; Evidence for reproductive toxicity in experimental animals. CAS-No.: 10287-53-3, 75980-60-8 Repr. 1B

### STOT-single exposure:

CAS-No.: 75980-60-8 NOAEL: 100 mg/kg bw/day (based on the components)

### STOT-repeated exposure:

CAS-No.: 72869-86-4 NOAEL: 100 mg/kg bw/day (based on the components)

### Aspiration hazard:

May cause respiratory irritation.

## SECTION 12: Ecological information

### 12.1. Toxicity

There are no data available on the preparation/mixture itself.

#### Ecotoxicity

| Substance:                                                                                     | CAS-No.:   | Ecotoxicity                                                                                                             |
|------------------------------------------------------------------------------------------------|------------|-------------------------------------------------------------------------------------------------------------------------|
| Diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide                                                | 75980-60-8 | EC50 (Daphnia, 48 h) 3,53 mg/L<br>EC50 (algae, 72 h) > 2,01 mg/L<br>EC50 (microorganisms, 3h) > 1000 mg/kg              |
| 3,6,9-trioxaundecamethylene dimethacrylate                                                     |            | EC50 (Daphnia, 48 h) 391 mg/L<br>EC50 (algae, 72 h) 68 mg/L                                                             |
| Ethyl 4-dimethylaminobenzoate                                                                  | 10287-53-3 | LC50 (fish, 96 h) 1,9 mg/L<br>NOEC (fish, 96h) 1,2 mg/L<br>LC50 (Daphnia, 48h) 4,5 mg/L                                 |
| 7,7,9(or 7,9,9)-trimethyl-4,13-dioxo-3,14-dioxa-5,12-diazahexadecane-1,16-diyl bismethacrylate | 72869-86-4 | LC50 (Daphnia) > 1,2 mg/L<br>EC50 (algae, 72 h) > 0,68 mg/l<br>LC50 (fish, 96 h) 10,1 mg/L<br>NOEC (algae, 72h) 21 mg/l |

### 12.2. Persistence and degradability

There are no data available on the preparation/mixture itself. Based on data of the components: CAS-No.: 75980-60-8, 72869-86-4, 10287-53-3 Not readily biodegradable (according to OECD criteria)

### 12.3. Bioaccumulative potential

There are no data available on the mixture itself. Log KOW: 1.39-3.39 (based on the components)

### 12.4. Mobility in soil

There are no data available on the preparation/mixture itself. Based on data of the components: CAS-No.: 10287-53-3, 72869-86-4 Koc: 560, 4516 @20 °C

### 12.5. Results of PBT and vPvB assessment

There are no data available on the mixture itself. The substances in the mixture do not meet the PBT/vPvB criteria according to REACH, annex XIII. (based on the components)

### 12.6 Endocrine disruptive effect



The substance/mixture contains no components  $\geq 0.1\%$  classified as endocrine disruptors under REACH Art. 57(f), or Reg. (EU) 2017/2100 or 2018/605.

## 12.7. Other adverse effects

There are no data available on the mixture itself.

## SECTION 13: Disposal considerations

### 13.1. Waste treatment methods

#### Appropriate disposal/Product:

Dispose of waste according to applicable legislation. Waste disposal according to directive 2008/98/EC, covering waste and dangerous waste.

Liquid, uncured resin: 08 01 11, 08 04 09;

Cured resin residues: 07 02 13, 20 01 39;

Cleaning waste: (Isopropanol, CL-1) 14 06 03; aqueous washing liquids 16 10 01.

#### Appropriate disposal / Package

Non-contaminated packages may be recycled.

Metal container: 15 01 04;

Plastic packaging: 15 01 02;

Special rules on packaging: 15 01 10 - Packaging that contains residues of hazardous substances or is contaminated by hazardous substances.

#### List of proposed waste codes / waste designations according to EWC / AVV

The allocation of waste identity numbers/waste descriptions must be carried out according to the EEC, specific to the industry and process.

## SECTION 14: Transport information

### 14.1. UN number

UN No.: -

### 14.2. UN proper shipping name

Land transport (ADR/RID)

-  
-

Sea transport (IMDG), Air transport (ICAO-TI / IATA-DGR)

-  
-

### 14.3. Transport hazard class(es)

Hazard label(s) / Label: -

Classification code: / Classification Code: -

### 14.4. Packing group

Packing group/ Packing Group: -

### 14.5. Environmental hazards

ADR/RID / IMDG / ICAO-TI / IATA-DGR:  
Marine pollutant:

Yes  
☐  
☐

No  
☒  
☒

#### 14.6. Special precautions for user

##### Land transport (ADR/RID)

|                     |   |                          |   |
|---------------------|---|--------------------------|---|
| transport category: | - | tunnel restriction code: | - |
| Special provisions: | - | Limited quantity (LQ):   | - |

##### Sea transport (IMDG)

EmS-No: -

Special provisions: - Limited quantity (LQ): -

#### 14.7. Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Remark No bulk transport by sea is intended.

### SECTION 15: Regulatory information

#### 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

##### EU legislation

Information on Regulation (EC) No 166/2006 establishing a European Pollutant Release and Transfer Register:

Substance / mixture / product / ingredients not listed

**Regulation (EC) No. 1005/2009 on substances that lead to the depletion of the ozone layer:**

REGULATION (EU) 2024/590: not applicable , Substance / mixture / product / ingredients not listed

**Regulation (EC) No. 648/2004 (Detergents regulation)**

not applicable

**Regulation (EC) No 850/2004 [POP-Regulation]:**

REGULATION (EU) 2019/102: Substance / mixture / product / ingredients not listed

**Regulation (EU) No 649/2012 on the export and import of dangerous chemicals:**

Substance / mixture / product / ingredients not listed

**Use restriction according to REACH annex XVII, no.:**

Use restriction according to REACH annex XVII, no.: 75 not applicable, 78 Nr.5

##### National regulations

Observe in addition any national regulations!

##### Restrictions of occupation

Observe restrictions to employment for juvenils according to the 'juvenile work protection guideline' (94/33/EC). Observe employment restrictions under the Maternity Protection Directive (92/85/EEC) for expectant or nursing mothers.

##### Other regulations, restrictions and prohibition regulations

This mixture contains the following substances of very high concern (SVHC) which are included in the Candidate List according to Article 59 of REACH: CAS No. 75980-60-8 Use restriction according to REACH annex XVII, no.: not applicable not applicable ; CA Proposition 65 : Substance(s) not listed. ; RoHS Directive 2011/65/EU : not relevant

#### 15.2. Chemical Safety Assessment

For this preparation a chemical safety assessment has been carried out. -

For this substance a chemical safety assessment has not been carried out.

## SECTION 16: Other information

### Relevant H- and EUH-phrases (Number and full text):

#### Hazard statements

H317 May cause an allergic skin reaction.  
H360 May damage fertility or the unborn child.

#### Training advice

Observe instructions for use.

#### Recommended restrictions of use:

Reserved for industrial and professional use. The product should only be handled by persons over the age of 18, who were informed sufficiently about the dangerous nature of the product and about the necessary safety precautions. Do not use for private purposes (household).

#### Further remarks:

The information provided in this Safety Data Sheet is based on our current knowledge at the time of publication. It is intended to offer guidance for the safe handling of the product mentioned herein during storage, processing, transport, and disposal. These details are not transferable to other products. If the product is blended, mixed, or processed with other materials, or subjected to any treatment, the information in this Safety Data Sheet may not apply to the resulting new material, unless explicitly stated otherwise.

#### Documentation of changes:

Changes from version 2.1:

Revision of all sections.

#### Key literature references and sources for data

Data arise from reference works and literature.

#### Abbreviations and acronyms

AC: Article category  
ACGIH: American Conference of Governmental Industrial Hygienists  
ADN: European agreement concerning the international carriage of dangerous goods by inland waterways (Accord européen relatif au transport international des marchandises dangereuses par voies de navigation intérieures)  
ADR: European agreement concerning the international carriage of dangerous goods by road (Accord européen relatif au transport international des marchandises dangereuses par route)  
AGW: Occupational exposure limit  
AOX: Adsorbable organically bound halogens  
Bw/KG: Body weight  
CAS: Chemical Abstracts Service  
CMR: Substances classified as carcinogenic, mutagenic or toxic for reproduction  
CSR: Chemical safety report  
DIN: German Institute for Standardization (Deutsches Institut für Normung)  
DNEL: Derived no-effect level  
DPD: Dangerous Preparations Directive (Richtlinie 1999/45/EC)  
DSD: Dangerous Substances Directive (Richtlinie 67/548/EWG)  
DU: Downstream user  
EC50: Effective concentration 50%  
ECHA: European Chemicals Agency  
EN: European standard  
EWC/EWL: European Waste Catalogue  
GHS: Globally harmonized system of classification and labelling of chemicals  
IATA: International Air Transport Association  
IBC: Intermediate bulk container  
ICAO: International Civil Aviation Organization  
IMDG Code: International Maritime Dangerous Goods Code  
IMO: International Maritime Organization  
ISO: International Organization for Standardization  
LC50: Lethal concentration 50%  
LD50: Lethal dose 50%  
LEV: Local exhaust ventilation  
LOAEL: Lowest observed adverse effect level  
LOEL: Lowest observed effect level  
MAK: Maximum workplace concentration (Maximale Arbeitsplatzkonzentration – DFG)  
n.a.: Not applicable  
n.b.: Not determined  
NOAEC: No observed adverse effect concentration  
NOAEL: No observed adverse effect level  
NOEC: No observed effect concentration  
OEL: Occupational exposure limit  
PBT: Persistent, bioaccumulative and toxic  
PNEC: Predicted no-effect concentration  
POP: Persistent organic pollutant  
PPE/PSA: Personal protective equipment  
REACH: Registration, Evaluation, Authorisation and Restriction of Chemicals  
RID: Regulations concerning the international carriage of dangerous goods by rail (Règlement international concernant le transport des marchandises dangereuses par chemin de fer)  
STEL: Short-term exposure limit  
SVHC: Substance of very high concern  
TLV: Threshold limit value  
UN: United Nations  
VOC: Volatile organic compounds  
vPvB: Very persistent and very bioaccumulative  
dw: Dry weight